


2015 -2016



DETERMINATION OF MELTING POINT

Assistant Lecturer **Sahar Mohammed**
Assistant Lecturer **Wid Kadhim**



Physical Properties of Organic Compounds

The physical properties of a compound include such things as its,

Color , Odor , Refractive index, Density, Solubility, Melting point & Boiling point.

The exact values of the physical properties of a cpd. depend on its **molecular structure**.

 in laboratory often depends on making a good prediction of physical properties from the cpd. structure.

The physical properties of a cpd. depend largely upon which kind of bonds hold it's atoms together in a molecule.

There are two kinds of bonds:

- 1- Ionic Bonds, formed by the transfer of electrons.
- 2- Covalent Bonds, formed by the sharing of electrons.


Phases of Matter

Matter ordinarily exists in four phases:

Solid , **Liquid** , **Gas** and **Plasma** .

Plasma is a high – temperature phase not encountered in the typical organic chemistry laboratory.


As the environment around the materials is changed, for example by varying the temp. or pressure, the systems undergo a phase transition, that is the change from one state of matter into another.



There are four Basic types of Intermolecular Forces :



1. Ion - dipole:

Ion is attracted to polar molecule via this force as ex. (NaCl in water).



2. Dipole – dipole:

Polar molecules attracted to each other as forces between HCl molecules.




3. Dipole - induce dipole:

Polar molecules attracted to nonpolar molecules as for ex.

(Oxygen in water)

by dipole – induce dipole interaction .



4. Induce dipole -induce dipole:


London dispersion forces also called van der Waal's forces
Nonpolar molecules attraction for each other due to electron distortion.

VAN DER WAALS' FORCES (VDW)
DIAGRAM

KEY

+ POSITIVE NUCLEUS

— NEGATIVE CHARGED ELECTRON CLOUD





When two atoms come within 5 nanometers of each other, there will be a slight interaction between them, thus causing polarity and a slight attraction.

Hydrogen bonding :

Hydrogen bonding is an especially strong kind of dipole – dipole attraction in which a hydrogen atom serves as a bridge between 2 electronegative atoms, holding one by a covalent bond & the other by purely electrostatic forces.

For H-bonding to be important, both electronegative atoms must come from the group: F, O, N.



Melting Point:

The melting point of a solid crystalline compound is the temperature at which the solid and the liquid phase of the cpd. are in equilibrium at a certain pressure usually 1 atmosphere.

Or,

It's the temp. at which the solid begins to change into liquid under a pressure of 1 atmosphere.

The melting point range : is the difference between the temp. at which the sample begins to melt (T_1) & the temp. at which the sample is completely melted (T_2) .

Melting point rang = $T_2 - T_1$

Organic chemist used melting point to:

1- Get an indication of the purity of cpd.

Sample	Pure	Impure	Decomposed
Melting point	Sharpe	Not sharpen	Not sharpen
Melting point range	0.5 – 1 °C	Wider range	Wider range

2- Help in identification of organic cpd.

Compound name	Compound formula	Compound Melting point
Sodium chloride	NaCl	801 °C
Sodium acetate	CH ₃ COO Na	324 °C
Glacial acetic acid	CH ₃ COOH	16.6 °C

Mixed Melting Point:

How can you know whether the 2 samples are the same or different?

Answer:

We mix them and measure the m.p. for the resultant mixture,
If,

m.p. of mixture < m.p. of original & wider m.p. range → different cpd.

m.p. of mixture = m.p. of original & same m.p. range → same cpd.




A
m.p. = 120 - 122 °C




B
m.p. = 120 - 122 °C

General Technique for m.p. Determination:

Capillary Tube method.



**Tube tilted to prevent
water condensation
inside**



**Stages for correct
sealing**



**Wrong !
Tube has been
overheated**



Filling a capillary tube:

The open end of the tube is pressed gently & repeatedly into the powdered sample several times.



The sample can be moved to the sealed end by turning the tube over & tapping it on the bench.





Any compound being loaded into a capillary tube must be:

1- Fully dried.

2- Homogeneous.

3- In powdered form.


Melting Point Apparatus:



(a)



(b)



(c)



(d)



(e)

a- Hot - oil m.p. bath.

c- Long necked flask.

e- Digital m.p. measuring apparatus.

b- Thiele apparatus.

d- Mel-temp apparatus.

Number of Experiment: **I**

Name of Experiment: **Determination of the
Melting Point
of an unknown sample.**

Aim of Experiment:


- 1- Identification of an unknown cpd. using it's m.p.
- 2- Determination of purity of a cpd. using the m.p.
as a physical property.

Procedure:

- 1- Preparation of the unknown sample as mentioned before.

2- The capillary tube is attached to a thermometer by a rubber ring in such a way that the closed end is attached to the bottom of the thermometer's bulb.

3- Then both of them are placed in an oil - bath , (the rubber ring should be above the surface of the oil - bath).




4- Heating is started gradually.

5- The range between the temperature at which the powdered solid inside the capillary tube begins to liquefy (T_1) & the temperature at which a clear liquid is observed inside the capillary tube (T_2) is recorded as the observed m.p. range.

$$\text{m.p. rang} = T_2 - T_1$$

$$\text{m.p.} = (T_1 - T_2)$$



Notes:

- * For oil bath, glycerin or paraffin can be used, but glycerin is preferred since it's safer than paraffin.
- * The temperature of the oil bath should rise very slowly.
- * Gentle heating with continuous stirring is necessary
- * After completing the experiment, the oil liquid must be allowed to cool before using it again.



Salt is added on snow covered roads?

Ice melts at ---- ?
What happens to the
freezing point if salt
is added?

Among these 3 samples of the same
cpd. , which of them is the purest one?



m.p.
115 - 119°C



m.p.
118 - 120 °C



m.p.
121- 122°C

Two of these bottles contain benzoic acid and one m-nitrophenylacetic acid.



1

m.p.=120 -122 °C



2

m.p.=120 -122°C



3

m.p.=120 -122°C

Run mixed melting points to confirm identification :

Mixed 1 and 2 m.p. = 120 - 122°C

Mixed 1 and 3 m.p. = 114 - 117°C

Mixed 2 and 3 m.p. = 115 - 118°C