

2015 - 2016



## **Back Titration**

# **Quantitative Assay of Aspirin Tablets**

*Assistant Lecturer*

**Sahar Mohammed Shakir**

*Assistant Lecturer*

**Abdul Hafeedh Hameed**

*Assistant Lecturer*

**Ali Basim**



# **Methods of Analytical Chemistry is of two types:**

## **1- Qualitative Analysis:**

**It determines the presence or absence of a particular compound, but not the mass or concentration.**

**By definition , qualitative analysis do not measure quantity.**

## **2- Quantitative Analysis:**

**It determines how much of each component , or of specified components is present in a given sample.**

### ***Methods of Quantitative Chemical Analysis :***

**1-Volumetric ( Titrimetric ) analysis.**

**2-Gravimetric analysis.**

**3-Spectrophotometric analysis.**

# **Requirements For a Titremetric Assay:**

**1.**

**The reaction can be represented by a chemical equation.**

**2.**

**The reaction should be relatively fast.**

**3.**

**The reaction should be complete & irreversible.**

**4.**

**The end point should be easily detected.**

## **Types of Titration:**

- 1- Forward titration (*direct titration*).
- 2- Back titration (*indirect titration*).

### **Back Titration:**

It includes the addition of an excess of a std. solution to a weighted amount of a sample and then the excess unreacted std. solution is determined by titration with another std. solution.


### **Back Titration Is Used For:**

- 1- Volatile substances, e.g.,  **$\text{NH}_3$** .
- 2- Insoluble or slightly soluble substances, e.g.  **$\text{CaCO}_3$**
- 3- Substances for which the quantitative reaction proceeds rapidly only in the presence of excess of reagent, e.g., **Lactic acid** & **Aspirin**.
- 4- Substances which decompose on heating, e.g. , **Formaldehyde**.

# Assay Of Aspirin


## Principle:

The determination of the amount of aspirin present in a tablet dosage form is done by alkaline hydrolysis of aspirin using  $N/2$  NaOH standard solution followed by back titrating of the excess unreacted alkali using  $N/2$  HCl std. solution & phenol red as indicator.



**Aspirin readily dissolved in dilute NaOH solution and hydrolyzed completely by heating for 10 minutes with an excess of a base.**


**Titration of the excess unreacted alkali with N/2 HCl std. solution using phenol red indicator**



**As in other quantitative determination involving boiling with a standard alkali, cooling and back titrating the excess, it's necessary to carry out a blank experiment without the aspirin**

**In order to:**

- 1- Minimize any error due to small unavoidable losses.**
- 2- Heating and cooling an alkaline liquid results in an apparent change in strength if certain indicators are used .**





**This change may be due to the interaction of the reagent with the glass or due to , the absorption of atmospheric CO<sub>2</sub> , CO<sub>2</sub> is rapidly absorbed by the hot alkaline solution to form sodium carbonate .**





**In the back titration with the standard acid the liberated CO<sub>2</sub> causes a color change of the indicator before the actual end point.**







# Phenol Red Indicator:


It's also known as phenolsulfonphthalein ( PSP ) is a pH indicator .





8.2  
Basic media  
Pink to red



6.4  
Acidic media  
Yellow





## Chemical Factor

1 mole of aspirin  $\equiv$  2 moles of NaOH

180 g of aspirin  $\equiv$  2 Eq. of NaOH

180/2 g of aspirin  $\equiv$  1000 ml of 1N NaOH

90/2 g of aspirin  $\equiv$  1000 ml of N/2 NaOH

0.045 g of aspirin  $\equiv$  1ml of N/2 NaOH

Volumes should be corrected to 0.5 N since procedure require that.



**A student had crushed 20 aspirin tablets ( 0.3 g aspirin each ) , recorded the weight of the powder as 6.5 g & he quantitatively analyzed aspirin by back titration .**

**Using 0.5 g of crushed powdered sample, N/2 NaOH & N/2 HCl std. sln. He recorded the results in the following data table:**

	<i>Sample exp.</i>	<i>Blank exp.</i>
<b>Weight of the powdered aspirin used</b>	<b>0.5 g</b>	<b>_____</b>
<b>Volume of 0.49 N NaOH used</b>	<b>30 ml</b>	<b>30 ml</b>
<b>Volume of 0.51 N HCL consumed</b>	<b>19 ml</b>	<b>28 ml</b>

**Each 0.045 g of aspirin  $\equiv$  1 ml of 0.5 N NaOH std. sln.  
 $\equiv$  1 ml of 0.5 N HCl std. sln.**

**1- Calculate the average weight of an individual aspirin tablet?**

**Answer:**

**The average wt. of an individual tablet =  $6.5 / 20$   
= 0.325 g**

**2- How many grams of aspirin ( pure ) is present in the student's sample?**

**Answer:**

**1<sup>st</sup> We should correct the normality of the used HCl to 0.5 N**

$$N_1 * V_1 = N_2 * V_2$$

$$0.51 * 19 = 0.5 * V_2$$

**$V_2 = 19.38$  ml of 0.5 N HCl for back titration  
( reacted with the excess NaOH std. solution )**

**2<sup>nd</sup> We should also correct the volume of blank**

$$N * V = N_3 * V_3$$

$$0.51 * 28 = 0.5 * V_3$$

**$V_3 = 28.56$  ml of 0.5 N HCl used for blank titration.**

**3<sup>rd</sup>** Calculate the volume of 0.5 N HCl reacted with pure aspirin indirectly.

$$V_3 - V_2 = 28.56 - 19.38 = 9.18 \text{ ml of } 0.5 \text{ N HCl solution.}$$

**4<sup>th</sup>** Calculate the amount of pure aspirin present in the sample by using the calculated **Chemical Factor**.

*Each 1ml of 0.5 N HCl is equivalent to 0.045 g of Aspirin.*

$$9.18 * 0.045 = 0.413 \text{ g of aspirin in the sample.}$$

# Discussion Questions

- 1- Why did you use your burette and not a graduated cylinder to add the excess NaOH standard solution?**
- 2- What is the definition of:**
  - a- Titration?**
  - b- Equivalence point?**
  - c- End point?**
  - d- Standard solution?**
  - e- Indicator?**
  - f- Molarity & Normality?**