

2015 - 2016



# Solvent Extraction

*Assistant Lec.*


**Sahar Mohammed**

*Assistant Lec.*

**Wid Kadhim**

# Extraction:

It's a method of **separation & purification** of organic compd.s depends on the ability of the compd. to dissolve in two immiscible solvents, e.g.,  $\text{H}_2\text{O}$  &  $\text{CH}_3\text{Cl}$ .



**Immiscible solvents  
are mixtures of liquids  
Insoluble in each other**

## **Extraction with immiscible solvents is generally employed for:**

- 1- Isolation of dissolved cpd.s from solution.**
- 2- Isolation of solid cpd.s from mixtures.**
- 3- Removal of undesirable impurities From mixtures ( *washing* ).**
- 4- Sometimes it's used in organic chemistry laboratories to remove an organic cpd . from a solution when the use of distillation is not convenient.**


## **Choosing a solvent for extraction:**

- 1- Readily dissolve the cpd. to be extracted.**
- 2- Have a low boiling point so it can be readily removed.**
- 3- Not react with solute or other solvent.**
- 4- Not be flammable or toxic.**
- 5- Show little or no water solubility ( immiscible with water ) .**
- 6- Be inexpensive.**

**No solvent meets all these criteria,  
for ex., ether,**


**is probably the most common solvent  
used for extraction but it is flammable**

**Diethyl ether**, Ether,  
 $\text{C}_4\text{H}_{10}\text{O}$ , Colorless liquid,  
Boiling point  $34.6\text{ }^\circ\text{C}$   
Density  $0.7134\text{ g / ml}$ ,  
Solubility in water  $69\text{g/L}$



Ether:

- 1-** Has a high solvating power for hydrocarbons & O containing cpds.
- 2-** Is highly volatile, b.p.  $34.6\text{ }^\circ\text{C}$ .  
So that it can be easily removed from the extract at low temp.  
Thus even highly sensitive cpds are not likely to decompose.
- 3-** Is very slightly soluble in water.






**Ether**, is used extensively as an extracting solvent. It's very slightly soluble in water & it's efficiency in use can be  $\uparrow$  by the addition of a small amount of an ionizable salt,  $\text{NaCl}$ , to the water layer. This lead to an  $\uparrow$  in the polarity of this soln. that result in a  $\downarrow$  in the solubility of a non polar cpd. This is known as "salting out " process.

### SALTING OUT,

In aqueous soln.s of organic molecules, salt is added to separate the organic material from the salty aqueous phase.



# Partition coefficient:

Partition coefficient is the ratio of concentrations of a compound in the two phases of a mixture of 2 immiscible liquids at equilibrium. Normally one of the solvents is aqueous while the 2nd is organic.



$$K = \frac{\text{Conc. of compound in organic solvent}}{\text{Conc. of compound in water}}$$

$$K = \frac{\text{Solubility of compound in organic solvent ( g/100ml )}}{\text{Solubility of compound in water ( g/100ml )}}$$

$$K = \frac{C_{org}}{C_w} = \frac{\text{Wt.}_{org} / V_{org}}{\text{Wt.}_w / V_w}$$

Where,

$C_{org}$  : Concentration of the solute in the organic layer.

$C_w$  : Concentration of the solute in the aqueous layer.

$\text{Wt.}_{org}$  : Weight of the solute in the organic layer.

$\text{Wt.}_w$  : Weight of the solute in the aqueous layer.

$V_{org}$  : Volume of the organic solvent.

$V_w$  : Volume of the aqueous solvent.



# **Name of experiment: Solvent Extraction**

## **Aim of experiment:**

**Separation of 4-hydroxybenzaldehyde from an unknown mixture containing NaCl and 4-hydroxybenzaldehyde.**

## **Properties of the cpd.s to be separated:**

### 4-Hydroxybenzaldehyde:

**Molecular formula:  $C_7H_6O_2$ , molar mass: 122.12 g/mol**


**Yellow-tan powder, Soluble in ether, Slightly soluble in water, Sublimes.**


### Sodium chloride:

**Molecular formula :  $NaCl$  ,**


**molar mass : 58.44 g/mol**

**Ionic cpd, Colorless crystals, Soluble in water.**






*Density of water = 1 g/ml*  
*Density of ether = 0.71 g/ml*



**1-** Transfer an unkn. sample to a separatory funnel; add 20 ml ether & 20 ml of D.W.



**3-** Leave for 5-10 min. to complete separation of the 2 layers. The stopper must be removed. Then separate the 2 layers.



**2-** Shake gently for 15 - 20 min. until no further pressure is released from the funnel stem.



**4-** Dry the ethereal layer

