

## Answer to problem on sodium hydroxide solution

- a) The volume of 1.11 N hydrochloric acid solution used in titrating sodium hydroxide sample was for total alkali regardless of whether sodium carbonate is present as impurity or not; *i. e.*, the same volume of hydrochloric acid solution would be consumed if sodium carbonate was precipitated as barium carbonate (explain this fact by the related balanced chemical equations).

$$\begin{array}{ccc} \text{HCl} & & \text{NaOH} \\ N_1 \times V_1 & = & N_2 \times V_2 \\ 1.11 \times 4.5 & = & N_2 \times 10 \\ N_2 & = & 0.4995 \text{ the normality of the prepared sodium hydroxide} \\ & & \text{solution} \end{array}$$

- b) Any solid sodium hydroxide sample taken from the original container would contain sodium carbonate as impurity by 3% w/w. Thus, you have to calculate how much sodium hydroxide was used to prepare the 10 milliliters- sample. This is calculated from the volume of hydrochloric acid used for total alkali:

$$\begin{array}{ccc} 1.11 \text{ N HCl} & & 1 \text{ N HCl} \\ N \times V & = & \hat{N} \times \hat{V} \\ 1.11 \times 4.5 & = & 1 \times \hat{V} \\ \hat{V} & = & 4.995 \text{ mL of 1 N hydrochloric acid solution would be required} \\ & & \text{for total alkali (V}_3\text{)} \end{array}$$

$$wt = V_3 \times ch. factor$$

$$wt = 4.995 \times 0.04$$

$wt = 0.1998$  g of total alkali (sodium hydroxide and the impurifying sodium carbonate) calculated as sodium hydroxide was used originally to prepare 10 mL of solution

$$\frac{3}{100} \times 0.1998 = 0.005994 \text{ g of sodium carbonate was present in the sample used to prepare 10 mL of solution}$$

$$wt = V_2 \times ch. factor$$

$$0.005994 = V_2 \times 0.053$$

$V_2 = 0.113$  mL of 1 *N* hydrochloric acid solution was consumed by the amount of sodium carbonate present in the sample assayed

$$\begin{array}{l} 1.11 \text{ N HCl} \quad 1 \text{ N HCl} \\ N \times V = \hat{N} \times \hat{V} \end{array}$$

$$1.11 \times V = 1 \times 0.113$$

$\hat{V} = 0.1$  mL of 1.11 *N* hydrochloric acid solution would be consumed by the amount of sodium carbonate present in the sample assayed

It is important to note that the percent w/w is related to solid samples. If you were asked to calculate the percent w/v, then this is related to the prepared solution as follows:

$$\text{for total alkali: } \frac{0.1998}{10} \times 100 = 1.998 \% w/v$$

$$\text{for sodium carbonate: } \frac{0.005994}{10} \times 100 = 0.05994 \% w/v$$

# Assay of Magnesium Sulfate

$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

Mwt.=246.5gm/mole


- ⦿ Contain not less than 99.5%  $\text{MgSO}_4$  calculated with reference of dried substance.
- ⦿ Description: colorless crystal, or white crystal powder.
- ⦿ Solubility: soluble at  $20^\circ\text{C}$  in 1.5 parts of water, sparingly soluble in alcohol 96% i.e. practically insoluble in 96% ethanol.
- ⦿ Uses: laxative agent.

○ Assay:

Each 1ml of M/20 EDTA  $\equiv$  0.012325gm  $MgSO_4 \cdot 7H_2O$

Note:

End point from *pink* to *blue* or *bluish violet*




- The technique involves titrating *metal ions* with a *complexing agent or chelating agent (Ligand)* and is commonly referred to as ***complexometric titration***. This method represents the analytical application of a complexation reaction.
- In this type of titration we use ***metal ion indicator***.
- **Chemical principle:**
- ❖ It is type of complexometric titration.




- ❖ Magnesium ions form relatively weak complex with (EDTA ;Ethylenediaminetetra-acetic acid).
- ❖ The stability of metal complex is pH dependent. Lower the pH of the solution, lesser would be the stability of complex (because more H<sup>+</sup> ions are available to compete with the metal ions for ligand)
- ❖ But at pH 10 the reaction is rapid and quantitative .

- ❖  $\text{NH}_3$ -  $\text{NH}_4\text{Cl}$  buffer is used to maintain this pH.
- ❖ EDTA which is polyamine carboxylic acid is one of the most frequently used chelating agent:



- ❖ EDTA form soluble, stable 1:1 chelates with metal ion.



❖ Because of the low water solubility of the free acid form of EDTA, the disodium salt (Disodium EDTA ) is used in preparation of std. sol.


⊙ **Indicator: Metal Indicator**

➤ Metal indicators are used to detect the end point in complexometric titration.

➤ They are dyes which produce one color in the presence of metal ions and a different color when they are absent.


➤ The color changes are the result of displacement of the H<sup>+</sup> ions by the metal.

➤ The metal indicator form 1:1 complex (MI) which must be weaker than the Metal-EDTA complex.





⊙ **Chemical factor:** it forms 1:1 chelate so;



0.012325gm  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O} \equiv 1\text{ ml of M/20 EDTA}$


⊙ **Chemical principle:**

- At the onset of titration, the buffered reaction media contains the MI complex and free M:




- When the titrant (EDTA) is added to the system, it will progressively chelate the free metal ions in the solution, and then a competitive reaction will take place, 1:1 metal ions and EDTA, since the MI complex is weaker than  $\text{M\_EDTA}$ .

- Finally, at the equivalent point, EDTA removes the last traces of metal ions from its colored complex with ind.(MI) to free metal ind.(I) of another color.



⦿ **Mordant black indicator:** it has blue color at pH10 and most of its complexes are reddish.

**Note:** In pH less than 6.3 and more than 11.5, the indicator and its complexes are reddish, so it is necessary to carry out the titration in the presence of **buffer at pH 10**.



○ Calculation:

$V \times \text{ch. factor} = \text{gm wt. of MgSO}_4 \cdot 7\text{H}_2\text{O}$   
 $\text{wt} \times 100/10 = \%w/v \text{ of unknown.}$

# Assay of Magnesium Sulfate

$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

Mwt.=246.5gm/mole


- ⦿ Contain not less than 99.5%  $\text{MgSO}_4$  calculated with reference of dried substance.
- ⦿ Description: colorless crystal, or white crystal powder.
- ⦿ Solubility: soluble at  $20^\circ\text{C}$  in 1.5 parts of water, sparingly soluble in alcohol 96% i.e. practically insoluble in 96% ethanol.
- ⦿ Uses: laxative agent.

○ Assay:

Each 1ml of M/20 EDTA  $\equiv$  0.012325gm  $MgSO_4 \cdot 7H_2O$

Note:

End point from *pink* to *blue* or *bluish violet*




- The technique involves titrating *metal ions* with a *complexing agent or chelating agent (Ligand)* and is commonly referred to as ***complexometric titration***. This method represents the analytical application of a complexation reaction.
- In this type of titration we use ***metal ion indicator***.
- **Chemical principle:**
- ❖ It is type of complexometric titration.




- ❖ Magnesium ions form relatively weak complex with (EDTA ;Ethylenediaminetetra-acetic acid).
- ❖ The stability of metal complex is pH dependent. Lower the pH of the solution, lesser would be the stability of complex (because more H<sup>+</sup> ions are available to compete with the metal ions for ligand)
- ❖ But at pH 10 the reaction is rapid and quantitative .



- ❖ NH<sub>3</sub>- NH<sub>4</sub>Cl buffer is used to maintain this pH.
- ❖ EDTA which is polyamine carboxylic acid is one of the most frequently used chelating agent:



- ❖ EDTA form soluble, stable 1:1 chelates with metal ion.



❖ Because of the low water solubility of the free acid form of EDTA, the disodium salt (Disodium EDTA ) is used in preparation of std. sol.


⊙ **Indicator: Metal Indicator**

➤ Metal indicators are used to detect the end point in complexometric titration.


➤ They are dyes which produce one color in the presence of metal ions and a different color when they are absent.

➤ The color changes are the result of displacement of the H<sup>+</sup> ions by the metal.

➤ The metal indicator form 1:1 complex (MI) which must be weaker than the Metal-EDTA complex.



○ **Chemical factor:** it forms 1:1 chelate so;



0.012325gm  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O} \equiv 1\text{ ml of M/20 EDTA}$


○ **Chemical principle:**

- At the onset of titration, the buffered reaction media contains the MI complex and free M:




- When the titrant (EDTA) is added to the system, it will progressively chelate the free metal ions in the solution, and then a competitive reaction will take place, 1:1 metal ions and EDTA, since the MI complex is weaker than  $\text{M\_EDTA}$ .

- Finally, at the equivalent point, EDTA removes the last traces of metal ions from its colored complex with ind.(MI) to free metal ind.(I) of another color.



⦿ **Mordant black indicator:** it has blue color at pH10 and most of its complexes are reddish.

**Note:** In pH less than 6.3 and more than 11.5, the indicator and its complexes are reddish, so it is necessary to carry out the titration in the presence of **buffer at pH 10**.



○ Calculation:

$V \times \text{ch. factor} = \text{gm wt. of MgSO}_4 \cdot 7\text{H}_2\text{O}$   
 $\text{wt} \times 100/10 = \%w/v \text{ of unknown.}$

# ASSAY OF MAGNESIUM HYDROXIDE IN ANTACID SUSPENSIONS (MAALOX)




- Maalox suspension contains two active ingredients, aluminium hydroxide and magnesium hydroxide. These are both medicines called antacids. They are used to reduce the acidity of the stomach's juices.
- Magnesium ion is complexed with EDTA. The solution is titrated with standard EDTA solution to the first permanent appearance of violet color.
- Mordant black is used as the indicator.






○ **Chemical principle:**



The overall reaction is given by :



## ○ Procedure:

note: We must first dilute the Maalox suspension before assay:

1. Transfer 10mL unknown to 250-mL conical flask.
2. Add 5mL of pH-10 buffer.
3. Add 10ml D.W.
4. Add about 0.1gm of *Mordant Black indicator*.
5. Titrate with M/40 EDTA solution until the color changes from pink to violet color.



- **Calculation:**
- Molecular Weight of  $\text{Mg(OH)}_2 = 58.3$
- $\text{mg of Mg(OH)}_2 \text{ Found} = M (\text{EDTA}) \times V (\text{EDTA})$   
in mL  $\times$  Mwt. of  $\text{Mg(OH)}_2$



- Complexometric titrations use reagents containing ligands that form complexes with metal ions. The indicator is a substance that forms a colored complex with the ions, but this complex must be weaker than that formed with the reagent. The metal ion complexes preferentially with the reagent, until at the end-point no metal-indicator complex remains, and its colour disappears.

An example is the determination of calcium and/or magnesium using EDTA (ethylenediaminetetra-acetic acid).





# ASSAY OF MAGNESIUM HYDROXIDE IN ANTACID SUSPENSIONS (MAALOX)



- Maalox suspension contains two active ingredients, aluminium hydroxide and magnesium hydroxide. These are both medicines called antacids. They are used to reduce the acidity of the stomach's juices.
- Magnesium ion is complexed with EDTA. The solution is titrated with standard EDTA solution to the first permanent appearance of violet color.
- Mordant black is used as the indicator.



○ **Chemical principle:**



The overall reaction is given by :



## ○ Procedure:

note: We must first dilute the Maalox suspension before assay:

1. Transfer 10mL unknown to 250-mL conical flask.
2. Add 5mL of pH-10 buffer.
3. Add 10ml D.W.
4. Add about 0.1gm of *Mordant Black indicator*.
5. Titrate with M/40 EDTA solution until the color changes from pink to violet color.





- **Calculation:**
- Molecular Weight of  $\text{Mg(OH)}_2 = 58.3$
- $\text{mg of Mg(OH)}_2 \text{ Found} = M (\text{EDTA}) \times V (\text{EDTA})$   
in mL  $\times$  Mwt. of  $\text{Mg(OH)}_2$



- **Complexometric titrations** use reagents containing ligands that form complexes with metal ions. The indicator is a substance that forms a colored complex with the ions, but this complex must be weaker than that formed with the reagent. The metal ion complexes preferentially with the reagent, until at the end-point no metal-indicator complex remains, and its colour disappears.

An example is the determination of calcium and/or magnesium using EDTA (ethylenediaminetetra-acetic acid).



# ASSAY OF BORAX

Borax or Sod. borate or  
Sod. Tetraborate

Chemical formula  $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$

Mwt. = 381.4 gm

- **Physical properties:** Borax is sod.tetraborate decahydrated, should be kept in well closed container, soluble in 20 parts of water.
- **Uses;** antiseptic and astringent in mouth paints.
- **Procedure:** 10mL of unknown

10mL of distilled water

2drops of methyl red indicator

Titrate with  $N/2$  HCl solution, the color of sol. changes from yellow  $\longrightarrow$  pink-orange (**1<sup>st</sup> end point**).

\* Record the volume of  $N/2$  HCl consumed.

Boil and cool the solution,

Add 4gm mannitol and 2drops of phenolphthalein ind. The color of sol. is pink


Now titrate with N/1 NaOH until the color of sol. changes from pink  $\rightarrow$  yellow  $\rightarrow$  pink(**2<sup>nd</sup> end point**).

\*Record the volume of N/1 NaOH consumed.

## 1<sup>st</sup> titration

N/2 HCl solution


10mL of unknown  
10mL of D.W  
2drops of methyl  
red indicator



## 2<sup>nd</sup> titration

N/1 NaOH solution

Boil , cool and add  
4gm mannitol and  
2drops of ph.ph ind.



- **Chemical principle:**

- ❖ Borax is a salt of weak acid with strong base.



- ❖ Boric acid ( $\text{H}_3\text{BO}_3$ ) will be liberated, it is very weak acid.

- ❖ Boil the solution to remove  $\text{CO}_2$ .


- ❖ Add mannitol because boric acid is weak acid to convert it into strong complex acid which is strong enough to be titrated directly with  $\text{NaOH}$ .

- ❖ Mannitol should be used in high conc. to prevent the hydrolysis of the complex.



sod. metaborate

❖ Mannitol should be used in high conc. to prevent the hydrolysis of this complex.



**mannityl boric acid**



- **Chemical factor:** *home work*
  - For borax with 0.5N HCl (0.09535 gm)
  - For borax with 1N NaOH (0.09535 gm)
  - For sod. carbonate with 0.5N HCl (0.0265gm)
  - For boric acid with 1N NaOH (0.06184gm)

Mwt. Boric acid=61.84 gm/mol

Mwt. Na<sub>2</sub>CO<sub>3</sub>= 106gm/mol

- **Calculations:**

**we must first correct the volume of HCl & NaOH**

1. If we have pure borax, the vol. of  $N/2$  HCl must be equal to the vol. of  $1N$  NaOH.

$V'_{\text{of } N/2 \text{ HCl}} = V'_{\text{of } 1N \text{ NaOH}}$  (the differences up to 0.3 mL)

take the average of the two volumes & calculate the % w/v of Borax.

2. Impure borax may contain  $\text{Na}_2\text{CO}_3$  or boric acid as impurities.

a) Mixture of **Borax &  $\text{Na}_2\text{CO}_3$**

$\text{Na}_2\text{CO}_3$  is basic; it can react with HCl, so both borax and carbonate will react with HCl



While NaOH react only with boric acid liberated from borax. Then after vol. correction;

$$V_1 \text{ of } N/2 \text{ HCl} > V_2 \text{ of } 1N \text{ NaOH}$$

$$V_1 - V_2 = \text{vol. of } N/2 \text{ HCl react with } \text{Na}_2\text{CO}_3$$

Calculate the % w/v for Borax and  $\text{Na}_2\text{CO}_3$ .

b) Mixture of **Borax & boric acid**



liberated boric acid



Liberated & free boric acid

So both the liberated and free boric acid react with NaOH, then

Correct the volumes and compare:

$V_2$  of 1N NaOH  $>$   $V_1$  of N/2 HCl

$V_2 - V_1 =$  vol. of 1N NaOH react with free boric acid

Calculate the % w/v for Borax and free boric acid.

# ASSAY OF BORAX

Borax or Sod. borate or  
Sod. Tetraborate

Chemical formula  $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$

Mwt. = 381.4 gm

- **Physical properties:** Borax is sod.tetraborate decahydrated, should be kept in well closed container, soluble in 20 parts of water.
- **Uses;** antiseptic and astringent in mouth paints.
- **Procedure:** 10mL of unknown

10mL of distilled water

2drops of methyl red indicator

Titrate with  $N/2$  HCl solution, the color of sol. changes from yellow  $\longrightarrow$  pink–orange(**1<sup>st</sup> end point**).

\*Record the volume of  $N/2$  HCl consumed.

Boil and cool the solution,

Add 4gm mannitol and 2drops of phenolphthalein ind. The color of sol. is pink


Now titrate with N/1 NaOH until the color of sol. changes from pink → yellow → pink(**2<sup>nd</sup> end point**).

\*Record the volume of N/1 NaOH consumed.

## 1<sup>st</sup> titration

N/2 HCl solution


10mL of unknown  
10mL of D.W  
2drops of methyl  
red indicator



## 2<sup>nd</sup> titration

N/1 NaOH solution

Boil , cool and add  
4gm mannitol and  
2drops of ph.ph ind.





- **Chemical principle:**

- ❖ Borax is a salt of weak acid with strong base.



- ❖ Boric acid ( $\text{H}_3\text{BO}_3$ ) will be liberated, it is very weak acid.

- ❖ Boil the solution to remove  $\text{CO}_2$ .


- ❖ Add mannitol because boric acid is weak acid to convert it into strong complex acid which is strong enough to be titrated directly with  $\text{NaOH}$ .

- ❖ Mannitol should be used in high conc. to prevent the hydrolysis of the complex.



sod. metaborate

❖ Mannitol should be used in high conc. to prevent the hydrolysis of this complex.



**mannityl boric acid**

- **Chemical factor:** *home work*
  - For borax with 0.5N HCl (0.09535 gm)
  - For borax with 1N NaOH (0.09535 gm)
  - For sod. carbonate with 0.5N HCl (0.0265gm)
  - For boric acid with 1N NaOH (0.06184gm)

Mwt. Boric acid=61.84 gm/mol

Mwt. Na<sub>2</sub>CO<sub>3</sub>= 106gm/mol

- **Calculations:**

**we must first correct the volume of HCl & NaOH**

1. If we have pure borax, the vol. of  $N/2$  HCl must be equal to the vol. of  $1N$  NaOH.

$V'_{\text{of } N/2 \text{ HCl}} = V'_{\text{of } 1N \text{ NaOH}}$  (the differences up to 0.3 mL)

take the average of the two volumes & calculate the % w/v of Borax.

2. Impure borax may contain  $\text{Na}_2\text{CO}_3$  or boric acid as impurities.

a) Mixture of **Borax &  $\text{Na}_2\text{CO}_3$**

$\text{Na}_2\text{CO}_3$  is basic; it can react with HCl, so both borax and carbonate will react with HCl



While NaOH react only with boric acid liberated from borax. Then after vol. correction;

$$V_1 \text{ of } N/2 \text{ HCl} > V_2 \text{ of } 1N \text{ NaOH}$$

$$V_1 - V_2 = \text{vol. of } N/2 \text{ HCl react with } \text{Na}_2\text{CO}_3$$

Calculate the % w/v for Borax and  $\text{Na}_2\text{CO}_3$ .

b) Mixture of **Borax & boric acid**



liberated boric acid



Liberated & free boric acid

So both the liberated and free boric acid react with NaOH, then

Correct the volumes and compare:

$V_2$  of 1N NaOH  $>$   $V_1$  of N/2 HCl

$V_2 - V_1 =$  vol. of 1N NaOH react with free boric acid

Calculate the % w/v for Borax and free boric acid.

# Assay of Chlorinated lime

Or bleaching powder

Or Ca. Chlorohypochlorite



# Assay of Chlorinated lime



- Contain not less than 30%w/w available chlorine.
- It is dull-white powder with characteristic odor, partially soluble in water and in 96% ethanol.
- Used as disinfectant.



## Procedure:

The flask used in the titration is *iodine flask*

$\approx 0.1\text{N Na}_2\text{S}_2\text{O}_3$  

The buret contains the titrant.

10ml of unknown  
1.5gm KI  
5ml of 33%v/v HAC


This flask contains the solution to be titrated and the indicator.



**note:** starch is added near the end point when the color of the solution changed from brown to yellow.

\* **End point is colorless**

- Chemical principle:



available chlorine

The *available chlorine* is the  $\text{Cl}_2$  which is liberated from the substance by the action of acid.

This liberated chlorine will displace an equivalent amount of  $\text{I}_2$  from  $\text{KI}$ .



The liberated I<sub>2</sub> is titrated against sod. thiosulfate using starch mucilage as indicator.



**Chemical factor:**

1ml of 0.1N Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>  $\equiv$  0.003545gm *available chlorine*.

**Calculation:**

Correct the sod. thiosulfate volume = V'

V' X ch. factor = gm of *available chlorine* in 10ml.

(Gm/10) X 100 = %w/v of *available chlorine*.

# Assay of Chlorinated lime

Or bleaching powder

Or Ca. Chlorohypochlorite



# Assay of Chlorinated lime



- Contain not less than 30%w/w available chlorine.
- It is dull-white powder with characteristic odor, partially soluble in water and in 96% ethanol.
- Used as disinfectant.

## Procedure:

The flask used in the titration is *iodine flask*

$\approx 0.1\text{N Na}_2\text{S}_2\text{O}_3$  

The buret contains the titrant.

10ml of unknown  
1.5gm KI  
5ml of 33%v/v HAC


This flask contains the solution to be titrated and the indicator.



**note:** starch is added near the end point when the color of the solution changed from brown to yellow.

\* **End point is colorless**

- Chemical principle:



available chlorine

The *available chlorine* is the  $\text{Cl}_2$  which is liberated from the substance by the action of acid.

This liberated chlorine will displace an equivalent amount of  $\text{I}_2$  from  $\text{KI}$ .



The liberated I<sub>2</sub> is titrated against sod. thiosulfate using starch mucilage as indicator.



**Chemical factor:**

1ml of 0.1N Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>  $\equiv$  0.003545gm *available chlorine*.

**Calculation:**

Correct the sod. thiosulfate volume = V'

V' X ch. factor = gm of *available chlorine* in 10ml.

(Gm/10) X 100 = %w/v of *available chlorine*.



# Assay of copper sulfate

Cupric sulfate  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$

***Mwt. = 249.68 gm***


- **Properties:** blue, fine prism or blue crystalline powder.
- **Solubility:** soluble at 20° c in 3parts of water and 3parts of glycerol, almost insoluble in 95% alcohol.  
It contains not less than 98.5% and not more than 101%w/w  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ .
- **Uses:** as fungicide.

- **Procedure:** the titration should be carried in stoppered flask.

During titration the color of the solution changed from **brown** to **faint yellow** then add 1ml starch and 1gm KSCN and complete titration.

$\approx 0.1N$  **Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>** →

10ml unknown  
30ml D.w.  
1.5 gm KI  
5ml acetic acid



- **Chemical principle:**

- ❖  $\text{CuSO}_4$  is treated with excess KI in acidic media.

oxidizing agent

acidic media



- ❖ The formed Cupric iodide ( $\text{CuI}_2$ ) is unstable.
- ❖  $\text{CuI}_2$  will break into cuprous iodide ( $\text{CuI}$  or  $\text{Cu}_2\text{I}_2$ ) and  $\text{I}_2$ .



unstable

yellow ppt.

- ❖ The liberated  $\text{I}_2$  is titrated with std.  $\text{Na}_2\text{S}_2\text{O}_3$  sol.



## • Notes:

1. KI is added to liberate I<sub>2</sub>.
2. Excess KI is added to dissolve I<sub>2</sub>.
3. HAC is added to make the medium slightly acidic.
4. In acidic media the liberated I<sub>2</sub> is adsorbed by the ppted CuI.
5. KSCN (pot. thiocyanate) is added near the end point.



this KI will return to  
react with CuSO<sub>4</sub> → → I<sub>2</sub>

- **Chemical factor:**

1ml of 0.1N  $\text{Na}_2\text{S}_2\text{O}_3 \equiv 0.02497\text{gm CuSO}_4 \cdot 5\text{H}_2\text{O}$

- **Calculation:**

- ❑ correct the volume of  $\text{Na}_2\text{S}_2\text{O}_3 = V'$
- ❑  $V' \times \text{Ch. Factor} = \text{gm wt. of CuSO}_4 \cdot 5\text{H}_2\text{O in 10ml}$
- ❑  $(\text{Gm}/10) \times 100 = \%w/v \text{ of CuSO}_4 \cdot 5\text{H}_2\text{O in the unknown}$

# Assay of copper sulfate

Cupric sulfate  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$

***Mwt. = 249.68 gm***

- **Properties:** blue, fine prism or blue crystalline powder.
- **Solubility:** soluble at 20° c in 3parts of water and 3parts of glycerol, almost insoluble in 95% alcohol.  
It contains not less than 98.5% and not more than 101%w/w  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ .
- **Uses:** as fungicide.




- **Procedure:** the titration should be carried in stoppered flask.

During titration the color of the solution changed from **brown** to **faint yellow** then add 1ml starch and 1gm KSCN and complete titration.

$\approx 0.1N$  **Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>**  $\longrightarrow$

10ml unknown  
30ml D.w.  
1.5 gm KI  
5ml acetic acid



- **Chemical principle:**

- ❖  $\text{CuSO}_4$  is treated with excess KI in acidic media.

oxidizing agent

acidic media



- ❖ The formed Cupric iodide ( $\text{CuI}_2$ ) is unstable.
- ❖  $\text{CuI}_2$  will break into cuprous iodide ( $\text{CuI}$  or  $\text{Cu}_2\text{I}_2$ ) and  $\text{I}_2$ .



unstable


yellow ppt.

- ❖ The liberated  $\text{I}_2$  is titrated with std.  $\text{Na}_2\text{S}_2\text{O}_3$  sol.



## • Notes:

1. KI is added to liberate I<sub>2</sub>.
2. Excess KI is added to dissolve I<sub>2</sub>.
3. HAC is added to make the medium slightly acidic.
4. In acidic media the liberated I<sub>2</sub> is adsorbed by the ppted CuI.
5. KSCN (pot. thiocyanate) is added near the end point.



this KI will return to  
react with CuSO<sub>4</sub> → → I<sub>2</sub>

- **Chemical factor:**

1ml of 0.1N  $\text{Na}_2\text{S}_2\text{O}_3 \equiv 0.02497\text{gm CuSO}_4 \cdot 5\text{H}_2\text{O}$

- **Calculation:**

- ❑ correct the volume of  $\text{Na}_2\text{S}_2\text{O}_3 = V'$
- ❑  $V' \times \text{Ch. Factor} = \text{gm wt. of CuSO}_4 \cdot 5\text{H}_2\text{O in 10ml}$
- ❑  $(\text{Gm}/10) \times 100 = \%w/v \text{ of CuSO}_4 \cdot 5\text{H}_2\text{O in the unknown}$

# Assay of Ferrous Sulfate



$M_{\text{wt.}} = 278 \text{ gm/mole}$


- **Appearance:** The blue-green heptahydrate is the most common form of ferrous sulfate ,it is bluish green crystals, freely soluble in water.
- **Content:** Ferrous Sulfate (crystal) contains 98.0—104.0% of ferrous sulfate heptahydrate ( $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ )
- **Used** in the prevention and treatment of iron deficiency.

- **Procedure:**

*H<sub>2</sub>SO<sub>4</sub> should be added to maintain a constant H ions concentration as explained*

$\approx 0.1\text{N KMnO}_4$  

*Unknown powder  
20ml of diluted H<sub>2</sub>SO<sub>4</sub>*




The buret contains the titrant.

This flask contains the solution to be titrated and the indicator.

*End point: colorless* 

- **Chemical principle:**



- Ferrous sulfate is a *reducing agent*.
- It will be titrated against standard solution of potassium permanganate (*oxidizing agent*) in acidic media.
- **Note:** you must write the oxidation and reduction equations by each.
- **Chemical factor:**

Derive the ch. Factor from the previous chemical equation.

$$1\text{ml. of } 0.1\text{N KMnO}_4 \equiv 0.0278\text{gm of FeSO}_4$$



- **Calculation:**

Correct the volume of  $\text{KMnO}_4$  consumed =  $V'$

Then multiply the corrected volume ( $V'$ ) by the Ch. Factor.

to get the weight;

$V' \times \text{ch. Factor} = \text{gm of FeSO}_4 \text{ in the unknown.}$

# Assay of Ferrous Sulfate



Mwt. = 278 gm/mole


- **Appearance:** The blue-green heptahydrate is the most common form of ferrous sulfate ,it is bluish green crystals, freely soluble in water.
- **Content:** Ferrous Sulfate (crystal) contains 98.0—104.0% of ferrous sulfate heptahydrate ( $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ )
- **Used** in the prevention and treatment of iron deficiency.

- **Procedure:**

*H<sub>2</sub>SO<sub>4</sub> should be added to maintain a constant H ions concentration as explained*

$\approx 0.1\text{N KMnO}_4$  

*Unknown powder*  
*20ml of diluted H<sub>2</sub>SO<sub>4</sub>*




The buret contains the titrant.

This flask contains the solution to be titrated and the indicator.

*End point: colorless* 

- **Chemical principle:**



- Ferrous sulfate is a *reducing agent*.
- It will be titrated against standard solution of potassium permanganate (*oxidizing agent*) in acidic media.
- **Note:** you must write the oxidation and reduction equations by each.
- **Chemical factor:**

Derive the ch. Factor from the previous chemical equation.

$$1\text{ml. of } 0.1\text{N KMnO}_4 \equiv 0.0278\text{gm of FeSO}_4$$

- **Calculation:**

Correct the volume of  $\text{KMnO}_4$  consumed =  $V'$

Then multiply the corrected volume ( $V'$ ) by the Ch. Factor.

to get the weight;

$V' \times \text{ch. Factor} = \text{gm of FeSO}_4 \text{ in the unknown.}$


# Assay of Hydrogen Peroxide



*H<sub>2</sub>O<sub>2</sub>: Mwt.=34.02 gm/mole*


*Different H<sub>2</sub>O<sub>2</sub> solution: 30%, 27%w/w and  
6%, 3%w/v*

- H<sub>2</sub>O<sub>2</sub> is slightly more viscous than water.
- Aqueous solution of H<sub>2</sub>O<sub>2</sub> colorless and odorless. Rapidly decomposes when contact with organic matter and certain metals; Cu, Fe, Mn and if allowed to stay in alkaline solution.






- In acidic solution  $\text{Fe}^{2+}$  is oxidized to  $\text{Fe}^{3+}$  ( $\text{H}_2\text{O}_2$  acting as an oxidizing agent):



- hydrogen peroxide acts as a reducing agent in its reaction with  $\text{KMnO}_4$  *in our experiment*.




- **Uses:**

It is a ***disinfectant*** or ***antiseptic*** for treating wounds for its antimicrobial activity.

it is an effective ***cleaning agent*** prevents infection transmission in the hospital environment.

# Procedure:

$\approx 0.1\text{N KMnO}_4$  

*10ml of unknown*  
*25ml of D.W.*  
*5ml of 50%v/v H<sub>2</sub>SO<sub>4</sub>* 



The buret contains the titrant.

This flask contains the solution to be titrated and the indicator.

**End point:** colorless  pink

- **Notes:**

1. We make the media acidic **to**;

- Prevent decomposition of  $\text{H}_2\text{O}_2$

- Prevent reaction of  $\text{KMnO}_4 \longrightarrow \text{MnO}_2$

- Increase oxidizing power of  $\text{KMnO}_4$

2. Titration should be done slowly


- **Chemical principle:**



*\*You must write the oxidation and reduction equations to calculate the eqwt. for both  $\text{KMnO}_4$  &  $\text{H}_2\text{O}_2$ .*

- **Chemical factor:**

According to the chemical equation;



*then complete the derivation of the Ch. Factor.*


**The Chemical factor is equal to:**



- **Calculation:**
- Correct the **volume** of  $\text{KMnO}_4$   
 $N \times V = N' \times V'$  (corrected)
- Multiply the corrected volume by the Ch. Factor.  
to get the **weight (wt.)** of  $\text{H}_2\text{O}_2$
- the **%w/v** of  $\text{H}_2\text{O}_2 = \text{wt.} / 10 \times 100$

- There is 2 methods to calculate the concentration of hydrogen peroxide in commercial products:

### 1. %W/V:

10ml of concentrated sol. was diluted to 200ml with D.W.  Then 20ml is taken from the dil. one and titrated with approximately 0.1N KMnO<sub>4</sub>

<u>conc.</u>	<u>diluted</u>	] x = $20 \times 10 / 200 = \underline{1\text{ml}}$ of conc. sol. (original sol.)
10	200	
x	20	

Correct the volume of KMnO<sub>4</sub> = **Vol.**

**Vol.** x ch. Factor = wt. of H<sub>2</sub>O<sub>2</sub> in 1ml.

**Then calculate %w/v**


## 2. Volume strength of the solution:

Decomposition occurs according to this equation:



Example: if the %w/v of  $\text{H}_2\text{O}_2$  is 8.5%w/v *i.e.* It contains 8.5 gm of  $\text{H}_2\text{O}_2$  in 100ml solution.

<u>wt.</u>	<u>Vol.</u>	
8.5	100ml	} X = 0.085 gm of $\text{H}_2\text{O}_2$ in 1ml
X	1ml (from the definition)	



So the volume strength of the sample is 28.

Assay of sodium  
hydroxide  
solution

NaOH solution



- From B.p, NaOH solution contains not less than 97.5% w/w of total alkali (as NaOH) and not more than 2.5 % w/w Na<sub>2</sub>CO<sub>3</sub>.
- **Assay:**
  - 10 ml of unknown(bulb pipette).
  - 25 ml of distilled water.
  - Add 2.5 ml of barium chloride solution.
  - Titrate with 1N HCl solution using 1-2 drops of phenolphthalein ind..
  - The first end point from pink → colorless(turbid)
  - To the turbid sol. add 5 drops of Bromophenol Blue ind. and complete titration with 1N HCl .
  - The second end point bluish violet → yellowish green

1N HCl solution



10 ml of unknown  
25ml of D.W  
2.5 ml of BaCl<sub>2</sub> solution  
1-2 drops of ph.ph. ind. (1<sup>st</sup> titration)  
5-6 drops of B.p.B ind. (2<sup>nd</sup> titration)



The buret contains  
the titrant.

This flask contains the  
solution to be titrated  
and the indicator.



- **1<sup>st</sup> titration:**



Why HCl do not react with BaCO<sub>3</sub>?

Why the end point is turbid?

- **2<sup>nd</sup> titration:**



***definition of chemical factor***: the weight of substance that is chemically equivalent to 1ml of std. solution.

- **Calculation of the *chemical factor*:**

a) From reaction of HCl with NaOH:

1Mwt of NaOH  $\equiv$  1 Mwt HCl

1 Mwt of NaOH  $\equiv$  1 eqwt HCL

1  $\times$  40 gm of NaOH  $\equiv$  1 liter of 1N HCl

40/1000 gm NaOH  $\equiv$  1ml of 1N HCl

0.04 gm of NaOH  $\equiv$  1ml of 1N HCl of total  
alkalinity calculated as NaOH(**chemical  
factor**)

b) From reaction of  $2\text{HCl}$  with  $\text{Na}_2\text{CO}_3$

$2\text{Mwt of HCl} \equiv 1\text{Mwt of BaCO}_3 \equiv 1\text{Mwt Na}_2\text{CO}_3$

$1\text{Mwt Na}_2\text{CO}_3 \equiv 2\text{Mwt of HCl}$

$1\text{Mwt Na}_2\text{CO}_3 \equiv 2 \text{ eqwt of HCl}$

$\frac{1}{2} \text{ Mwt Na}_2\text{CO}_3 \equiv 1 \text{ eqwt of HCl}$

$106/2 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{liter of } 1N \text{ HCL}$

$53 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{liter of } 1N \text{ HCL}$

$53/1000 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{ml of } 1N \text{ HCL}$

$0.053 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{ml of } 1N \text{ HCL (chemical factor)}$

- **Calculations :**

$V_1$  is the of HCl consumed in the 1<sup>st</sup> titration

$V_2$  is the of HCl consumed in the 2<sup>nd</sup> titration

$V_1 + V_2 = V_3$  total HCl consumed.

correct the  $V_3$  according to this equation:

$$V \times N = V' \times N'$$

Corrected  $V_3 \times 0.04 =$  gm wt. of total alkali

Then the % w/v of total alkali

Corrected  $V_2 \times 0.053 =$  gm wt. of  $\text{Na}_2\text{CO}_3$

Then the % w/v of  $\text{Na}_2\text{CO}_3$  in the unknown

**Preparation & Standardization of  $\approx 0.1$  N  
Potassium Permanganate ( $\text{KMnO}_4$ ) solution**

**Redox Reaction**



◦ Preparation of  $\approx 0.1N$   $KMnO_4$  solution:

- Dissolve 3.3 gm of solid  $KMnO_4$  in 1000ml. distilled water, in a flask.
- Heat the solution on a hotplate for 15min. at 70-80 °c.
- Cool the solution to room temp. and allow stand for 2 days in a dark place.
- Filter the solution through the glass wool or asbestos before use.

◦ **Standardization of  $\text{KMnO}_4$  solution:**

- ***Filter the prepared soln. through glass wool as said before.***
- ***Rinse and Fill the burette with filtered  $\text{KMnO}_4$  soln.***
- ***Weigh accurately about 200mg sod. oxalate ( $\text{Na}_2\text{C}_2\text{O}_4$ ), previously dried at  $110^\circ\text{C}$  to constant weight.***
- ***Dissolve it with 100ml of D.W. in a conical flask. Add 7ml. of conc.  $\text{H}_2\text{SO}_4$ .***
- ***Heat the soln. to about  $70^\circ\text{C}$ .***
- ***Titrate against  $\text{KMnO}_4$  soln. until a permanent pale pink color appear (persist for 15 sec.).***
- ***Record the volume of  $\text{KMnO}_4$  consumed.***

- **OXIDATION- REDUCTION REACTIONS:**

- ❖ Methods of analysis involve a change in valence of the reacting substances, thus there is always a transfer of electrons.
- ❖ Electrons lost by reducing agent are gained by oxidizing agent. The oxidizing agent will be reduced; gain of electrons and then decreasing in oxidation number, but the reducing agent will be oxidized; losing of electrons and then increasing in oxidation number.

- We can calculate the eqwt. of the reactants:

Mwt.

$$\text{Eq. wt} = \frac{\text{Mwt.}}{\text{no. of e- gained or lost by each molecule}}$$

no. of e- gained or lost by each molecule

- **Permanganate ion** ( $\text{MnO}_4^-$ ) acts as an oxidizing agent in acidic media, being reduced to  $\text{Mn}^{++}$ :



- Oxalate ion ( $\text{C}_2\text{O}_4^{=}$ ) acts as reducing agent:



- $\text{KMnO}_4$  is widely used in volumetric analysis as a standard and oxidizing agent because:
  1. It is strong oxidizing agent.
  2. The intense color of  $\text{KMnO}_4$  solution is sufficient to signal the end point in most titrations thus it is used as self indicator.
- On the other hand,
  - 1)  $\text{KMnO}_4$  is a strong oxidizing agent, then it will oxidize any organic matter present in D.W., so heating for 15min. is to accelerate the oxidation of these organic matter.

- 2)  $\text{KMnO}_4$  soln. is allowed to stand for 2 days to ensure the completion of the decomposition reaction.
- 3) It is filtered through asbestos to remove all traces of manganese dioxide ( $\text{MnO}_2$ ) .
- 4)  $\text{KMnO}_4$  is unstable in the presence of direct sunlight, organic matter ex. Filter paper.
- 5)  $\text{H}_2\text{SO}_4$  is added since  $\text{KMnO}_4$  is better oxidizing agent in acidic media,  $\text{H}_2\text{SO}_4$  is added to keep  $[\text{H}^+]$  ion conc. constant through out the titration.
- 6) Heat to about  $70^\circ\text{C}$  because the oxidation of sod. oxalate is rapid enough if the temp. above  $60^\circ\text{C}$  .

- Sod. oxalate, not oxalic acid, is used because:
  1. The salt can be obtained in very pure condition.
  2. Free from water of crystallization.
  3. It can be dried to 130 °c without decomposition.
  4. Stable to air.


The chemical equation for the std. of pot. Permanganate:



# Calculations

$$\frac{\text{KMnO}_4}{N \times V} = \frac{\text{Na}_2\text{C}_2\text{O}_4}{\frac{\text{WT}}{\text{eqwt}} \times 1000}$$

67gm/mole





# The preparation and standardization of 0.1 N sodium thiosulfate



## Preparation:

Dissolve about 26gm sod. thiosulfate and 200mg sod. carbonate in 1000ml of recently boiled and cooled D.W.

## Standardization:

Put into the stoppered flask the followings;

0.1–0.2 gm of potassium dichromate ( $K_2Cr_2O_7$ )

100ml of D.W.

3gm potassium iodide (KI)

2gm sod. Bicarbonate(  $NaHCO_3$  )


5ml conc. HCl

- ▶ Close the flask with the stopper and swirl to mix.
- ▶ Then keep the flask in dark place for 10min.

- ▶ Rinse the stopper and inner walls of the flask with D.W.
- ▶ Titrate the liberated  $I_2$  with sod. thiosulfate the color change from brown to faint yellow.
- ▶ Add 1 ml of starch indicator the color of the solution change to dark blue.
- ▶ Complete titration until the color change to green-blue.


≈ 0.1 N  $Na_2S_2O_3$

End point: dark blue → green-blue



## Notes for preparation:

- ▶ *Why do we use recently boiled and cooled water?*
  - a. to destroy sulfur bacteria
  - b. to expel CO<sub>2</sub>



- ▶ *Na<sub>2</sub>CO<sub>3</sub> is added to the solution as preservative*

## Chemical principle of standardization:

- ▶ Pot. Dichromate is an oxidizing agent



- ▶ The liberated  $\text{I}_2$  is titrated against  $\text{Na}_2\text{S}_2\text{O}_3$



sod. tetrathionate

- ▶ This indirect procedure is known as *iodometry*.

## Notes for standardization:


1. We use stoppered flask to prevent the loss of  $I_2$  .
2. KI is added to liberate  $I_2$ .
3. Excess KI is added, Why?



4. The standardization should be carried in acidic media and not in alkaline media:  
because alkaline media  $I_2$  will react with  $OH^-$  to give:



These ions will oxidize the thiosulfate to:



5. On the other hand in higher acidic media the excess KI will react with HCl to give HI, HI will be oxidized by atmospheric O<sub>2</sub> to I<sub>2</sub>:



excess

atmospheric O<sub>2</sub>



Sod. bicarbonate is added to react with formed HI:



6. We put the flask in dark place for 10min.?
7. Starch ind. is added at the end (when the color of the solution change to yellowish green) or we can said the ind. is added when the iodine concentration is low because  $\beta$ -amylase which the soluble form of starch react :  
$$\beta\text{-amylase} + \text{I}_2 + \text{I}^- \longrightarrow \text{blue colour complex}$$

the color change is reversible, the color being discharged (change to colorless) when iodine is reduced by  $\text{Na}_2\text{S}_2\text{O}_3$ .

This reversibility is reduced when the iodine conc. is high for this reason starch ind. should be added until most of  $\text{I}_2$  has been reduced by  $\text{Na}_2\text{S}_2\text{O}_3$  .



8. The end point is a pale green and not colorless why?

### Calculation:

$$\begin{array}{l} \text{Na}_2\text{S}_2\text{O}_3 \\ N \times V \\ \text{eqwt} = 49 \end{array} = \begin{array}{l} \text{K}_2\text{Cr}_2\text{O}_7 \\ (\text{Wt}/\text{eqwt}) \times 1000 \end{array}$$

# The preparation and standardization of 0.1 N sodium thiosulfate



## Preparation:

Dissolve about 26gm sod. thiosulfate and 200mg sod. carbonate in 1000ml of recently boiled and cooled D.W.

## Standardization:

Put into the stoppered flask the followings;

0.1–0.2 gm of potassium dichromate ( $K_2Cr_2O_7$ )

100ml of D.W.

3gm potassium iodide (KI)

2gm sod. Bicarbonate(  $NaHCO_3$  )


5ml conc. HCl

- ▶ Close the flask with the stopper and swirl to mix.
- ▶ Then keep the flask in dark place for 10min.

- ▶ Rinse the stopper and inner walls of the flask with D.W.
- ▶ Titrate the liberated  $I_2$  with sod. thiosulfate the color change from brown to faint yellow.
- ▶ Add 1 ml of starch indicator the color of the solution change to dark blue.
- ▶ Complete titration until the color change to green-blue.


≈ 0.1 N  $Na_2S_2O_3$

End point: dark blue → green-blue



## Notes for preparation:

- ▶ *Why do we use recently boiled and cooled water?*
  - a. to destroy sulfur bacteria
  - b. to expel CO<sub>2</sub>



- ▶ *Na<sub>2</sub>CO<sub>3</sub> is added to the solution as preservative*

## Chemical principle of standardization:

- ▶ Pot. Dichromate is an oxidizing agent



- ▶ The liberated  $\text{I}_2$  is titrated against  $\text{Na}_2\text{S}_2\text{O}_3$



sod. tetrathionate

- ▶ This indirect procedure is known as *iodometry*.

## Notes for standardization:


1. We use stoppered flask to prevent the loss of  $I_2$  .
2. KI is added to liberate  $I_2$ .
3. Excess KI is added, Why?



4. The standardization should be carried in acidic media and not in alkaline media:  
because alkaline media  $I_2$  will react with  $OH^-$  to give:



These ions will oxidize the thiosulfate to:



5. On the other hand in higher acidic media the excess KI will react with HCl to give HI, HI will be oxidized by atmospheric O<sub>2</sub> to I<sub>2</sub>:



excess

atmospheric O<sub>2</sub>



Sod. bicarbonate is added to react with formed HI:





6. We put the flask in dark place for 10min.?
7. Starch ind. is added at the end (when the color of the solution change to yellowish green) or we can said the ind. is added when the iodine concentration is low because  $\beta$ -amylase which the soluble form of starch react :  
$$\beta\text{-amylase} + \text{I}_2 + \text{I}^- \longrightarrow \text{blue colour complex}$$

the color change is reversible, the color being discharged (change to colorless) when iodine is reduced by  $\text{Na}_2\text{S}_2\text{O}_3$ .

This reversibility is reduced when the iodine conc. is high for this reason starch ind. should be added until most of  $\text{I}_2$  has been reduced by  $\text{Na}_2\text{S}_2\text{O}_3$  .

8. The end point is a pale green and not colorless why?

### Calculation:

$$\begin{array}{l} \text{Na}_2\text{S}_2\text{O}_3 \\ N \times V \\ \text{eqwt} = 49 \end{array} = \begin{array}{l} \text{K}_2\text{Cr}_2\text{O}_7 \\ (\text{Wt}/\text{eqwt}) \times 1000 \end{array}$$

# The preparation and assay of Aqueous Iodine Solution (Lugol's Solution)

- ▶ It contains 5% w/v *Iodine* ( limits 4.75– 5.25% w/v) and 10% w/v *Pot. Iodide* ( limits 9.5– 10.5% w/v).
- ▶ Used in pre-operative treatment of thyrotoxicosis.

- ▶ Preparation:

I<sub>2</sub> 50gm

KI 100gm

purified water sufficient to 1000ml

- \*Dissolve KI in a little amount of water then add the I<sub>2</sub> and shake well until the I<sub>2</sub> is completely dissolved, then add sufficient water to 1000ml using volumetric flask.

▶ Assay: Dilute 25ml of the original sol. to 100ml with distilled water, then: `

a) Assay of I<sub>2</sub>,


note: we must use stoppered flask

***Each 1 ml of N/10 Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> ≡ 0.01269 gm I<sub>2</sub>***

≈ N/10 Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> →

dark brown, titrate → yellow then add 1 ml starch  
the color change to dark blue complete tit.  
until the end point colorless

10ml of unknown  
10ml D.W




## b) Assay of KI,

\* titrate with  $M/20$   $KIO_3$  (Pot. Iodate) until **dark brown** sol. change to **pale brown**, add 1 ml of amaranth sol.; the color of sol. Change to **red** and then continue tit. Until the color of sol. change to **yellow -orange**

$M/20$   $KIO_3$  →

5 ml of unknown  
10 ml D.W  
20 ml conc. HCl



- ▶ From the no. of ml. of M/20  $\text{KIO}_3$  required subtract  $\frac{1}{4}$  no. of ml of N/10  $\text{Na}_2\text{S}_2\text{O}_3$  required in the assay of  $\text{I}_2$   
*each 1 ml of the remainder is  $\equiv$  0.0166 gm of KI*

▶ Chemical principle for I<sub>2</sub> assay:

I<sub>2</sub> will oxidize chemically equivalent amount of sod. thiosulfate.



▶ Chemical principle for KI assay:

KIO<sub>3</sub> (pot. Iodate) is fairly strong oxidizing agent, it can quantitatively react with both iodide I<sup>-</sup> and iodine I<sub>2</sub>.

\*If the conc. of HCl not exceed 1 N the reaction:





\*In the presence of conc. HCl exceeding 4N, I<sub>2</sub> produced by the previous reaction is further oxidized to iodine monochloride (ICl):



\* so the overall reaction between KIO<sub>3</sub> & KI in the presence of high conc. HCl >4N can be expressed by:



The presence of high conc. HCl is required to decrease :



stable cpd

- ▶ In this experiment we use amaranth indicator, it is a dye, dyes which are destroyed or changed in the presence of oxidizing titrant can be use as indicator.
- ▶ We can not use starch as indicator because starch\_ iodide complex can not be formed in high acid conc. used in the iodate titration.

▶ ***Chemical factor. Homework***

Each 1 ml of 0.1 N  $\text{Na}_2\text{S}_2\text{O}_3 \equiv 0.01269\text{gm}$  of Iodine

Each 1 ml of M/20  $\text{KIO}_3 \equiv 0.0166\text{gm}$  of KI

Each 1 ml of M/20  $\text{KIO}_3 \equiv 0.02538\text{gm}$  of iodine

▶ **Calculation:**


▶ **For I<sub>2</sub>:**

<u>original sol.</u>	<u>Vol.</u>	}	<i>X = 2.5ml of original sol.</i>
25	100		
X	10		


*V'* x ch.fact.  $0.01269 \times 100/2.5 = \%w/v$  of iodine in orig. sol.  
*V'* the corrected volume of sod. thiosulfate.

▶ **For KI:**

<u>original sol.</u>	<u>Vol.</u>	}	<i>X = 1.25ml of original sol.</i>
25	100		
X	5		



*Vol. of KIO<sub>3</sub> (req. for KI) x 0.0166 x 100/1.25 = %w/v of KI*

- ▶ It is preferred to use iodate ( $\text{KIO}_3$ ) of known M rather than of known N, it depends on the nature of reaction.
  - ▶  $\text{KIO}_3$  sol. used in the titration of sol. containing both iodide and iodine.
  - ▶ The equivalent of  $\text{KIO}_3$  in its reaction with KI differs from its equivalent when react with iodine.
- 

# The preparation and assay of Aqueous Iodine Solution (Lugol's Solution)

- ▶ It contains 5% w/v *Iodine* ( limits 4.75– 5.25% w/v) and 10% w/v *Pot. Iodide* ( limits 9.5– 10.5% w/v).
- ▶ Used in pre-operative treatment of thyrotoxicosis.

- ▶ Preparation:

I<sub>2</sub> 50gm

KI 100gm

purified water sufficient to 1000ml

- \*Dissolve KI in a little amount of water then add the I<sub>2</sub> and shake well until the I<sub>2</sub> is completely dissolved, then add sufficient water to 1000ml using volumetric flask.

▶ Assay: Dilute 25ml of the original sol. to 100ml with distilled water, then: `

a) Assay of I<sub>2</sub>,


note: we must use stoppered flask

***Each 1 ml of N/10 Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> ≡ 0.01269 gm I<sub>2</sub>***

≈ N/10 Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> →

dark brown, titrate → yellow then add 1 ml starch  
the color change to dark blue complete tit.  
until the end point colorless

10ml of unknown  
10ml D.W




## b) Assay of KI,

\* titrate with  $M/20$   $KIO_3$  (Pot. Iodate) until **dark brown sol. change to pale brown**, add 1 ml of amaranth sol.; the color of sol. Change to **red** and then continue tit. Until the color of sol. change to **yellow -orange**

$M/20$   $KIO_3$  →

5ml of unknown  
10ml D.W  
20ml conc. HCl





- ▶ From the no. of ml. of M/20  $\text{KIO}_3$  required subtract  $\frac{1}{4}$  no. of ml of N/10  $\text{Na}_2\text{S}_2\text{O}_3$  required in the assay of  $\text{I}_2$   
*each 1 ml of the remainder is  $\equiv$  0.0166 gm of KI*

▶ Chemical principle for I<sub>2</sub> assay:

I<sub>2</sub> will oxidize chemically equivalent amount of sod. thiosulfate.



▶ Chemical principle for KI assay:

KIO<sub>3</sub> (pot. Iodate) is fairly strong oxidizing agent, it can quantitatively react with both iodide I<sup>-</sup> and iodine I<sub>2</sub>.

\*If the conc. of HCl not exceed 1 N the reaction:



\*In the presence of conc. HCl exceeding 4N, I<sub>2</sub> produced by the previous reaction is further oxidized to iodine monochloride (ICl):



\* so the overall reaction between KIO<sub>3</sub> & KI in the presence of high conc. HCl >4N can be expressed by:



The presence of high conc. HCl is required to decrease :



stable cpd

- ▶ In this experiment we use amaranth indicator, it is a dye, dyes which are destroyed or changed in the presence of oxidizing titrant can be use as indicator.
- ▶ We can not use starch as indicator because starch\_ iodide complex can not be formed in high acid conc. used in the iodate titration.

▶ ***Chemical factor. Homework***

Each 1 ml of 0.1 N  $\text{Na}_2\text{S}_2\text{O}_3 \equiv 0.01269\text{gm}$  of Iodine

Each 1 ml of M/20  $\text{KIO}_3 \equiv 0.0166\text{gm}$  of KI

Each 1 ml of M/20  $\text{KIO}_3 \equiv 0.02538\text{gm}$  of iodine

▶ **Calculation:**

▶ **For I<sub>2</sub>:**

<u>original sol.</u>	<u>Vol.</u>	}	<i>X = 2.5ml of original sol.</i>
25	100		
X	10		

*V'* x ch.fact.  $0.01269 \times 100/2.5 = \%w/v$  of iodine in orig. sol.  
*V'* the corrected volume of sod. thiosulfate.

▶ **For KI:**

<u>original sol.</u>	<u>Vol.</u>	}	<i>X = 1.25ml of original sol.</i>
25	100		
X	5		

$$\begin{array}{l}
 \downarrow \\
 V_{\text{KIO}_3 \text{ M/20}} \text{ (req. for KI)} = V_{\text{KIO}_3 \text{ M/20}} \text{ (used in the total assay)} - \frac{1}{4} V'_{\text{Na}_2\text{S}_2\text{O}_3} \text{ 0.1N} \\
 \\
 \text{Vol. of KIO}_3 \text{ (req. for KI)} \times 0.0166 \times 100/1.25 = \%w/v \text{ of KI}
 \end{array}$$

- ▶ It is preferred to use iodate ( $\text{KIO}_3$ ) of known M rather than of known N, it depends on the nature of reaction.
  - ▶  $\text{KIO}_3$  sol. used in the titration of sol. containing both iodide and iodine.
  - ▶ The equivalent of  $\text{KIO}_3$  in its reaction with KI differs from its equivalent when react with iodine.
- 