

Assay of Magnesium Sulfate

$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

Mwt.=246.5gm/mole


- ⦿ Contain not less than 99.5% MgSO_4 calculated with reference of dried substance.
- ⦿ Description: colorless crystal, or white crystal powder.
- ⦿ Solubility: soluble at 20°C in 1.5 parts of water, sparingly soluble in alcohol 96% i.e. practically insoluble in 96% ethanol.
- ⦿ Uses: laxative agent.

○ Assay:

Each 1ml of M/20 EDTA \equiv 0.012325gm $MgSO_4 \cdot 7H_2O$

Note:

End point from *pink* to *blue* or *bluish violet*



- The technique involves titrating *metal ions* with a *complexing agent or chelating agent (Ligand)* and is commonly referred to as ***complexometric titration***. This method represents the analytical application of a complexation reaction.
- In this type of titration we use ***metal ion indicator***.
- **Chemical principle:**
- ❖ It is type of complexometric titration.



- ❖ Magnesium ions form relatively weak complex with (EDTA ;Ethylenediaminetetra-acetic acid).
- ❖ The stability of metal complex is pH dependent. Lower the pH of the solution, lesser would be the stability of complex (because more H⁺ ions are available to compete with the metal ions for ligand)
- ❖ But at pH 10 the reaction is rapid and quantitative .

❖ Because of the low water solubility of the free acid form of EDTA, the disodium salt (Disodium EDTA) is used in preparation of std. sol.


⊙ **Indicator: Metal Indicator**

➤ Metal indicators are used to detect the end point in complexometric titration.


➤ They are dyes which produce one color in the presence of metal ions and a different color when they are absent.

➤ The color changes are the result of displacement of the H⁺ ions by the metal.

➤ The metal indicator form 1:1 complex (MI) which must be weaker than the Metal-EDTA complex.



○ **Chemical factor:** it forms 1:1 chelate so;



0.012325gm $\text{MgSO}_4 \cdot 7\text{H}_2\text{O} \equiv 1\text{ ml of M/20 EDTA}$


○ **Chemical principle:**

- At the onset of titration, the buffered reaction media contains the MI complex and free M:




- When the titrant (EDTA) is added to the system, it will progressively chelate the free metal ions in the solution, and then a competitive reaction will take place, 1:1 metal ions and EDTA, since the MI complex is weaker than M_EDTA .

- Finally, at the equivalent point, EDTA removes the last traces of metal ions from its colored complex with ind.(MI) to free metal ind.(I) of another color.



- ⦿ **Mordant black indicator:** it has blue color at pH 10 and most of its complexes are reddish.

Note: In pH less than 6.3 and more than 11.5, the indicator and its complexes are reddish, so it is necessary to carry out the titration in the presence of **buffer at pH 10**.



○ Calculation:

$V \times \text{ch. factor} = \text{gm wt. of MgSO}_4 \cdot 7\text{H}_2\text{O}$
 $\text{wt} \times 100/10 = \%w/v \text{ of unknown.}$

ASSAY OF BORAX

Borax or Sod. borate or
Sod. Tetraborate

Chemical formula $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$

Mwt. = 381.4 gm

- **Physical properties:** Borax is sod.tetraborate decahydrated, should be kept in well closed container, soluble in 20 parts of water.
- **Uses;** antiseptic and astringent in mouth paints.
- **Procedure:** 10mL of unknown

10mL of distilled water

2drops of methyl red indicator

Titrate with $N/2$ HCl solution, the color of sol. changes from yellow \longrightarrow pink-orange (**1st end point**).

* Record the volume of $N/2$ HCl consumed.

Boil and cool the solution,

Add 4gm mannitol and 2drops of phenolphthalein ind. The color of sol. is pink


Now titrate with N/1 NaOH until the color of sol. changes from pink → yellow → pink(**2nd end point**).

*Record the volume of N/1 NaOH consumed.

1st titration

N/2 HCl solution


10mL of unknown
10mL of D.W
2drops of methyl
red indicator



2nd titration

N/1 NaOH solution

Boil , cool and add
4gm mannitol and
2drops of ph.ph ind.



- **Chemical principle:**

- ❖ Borax is a salt of weak acid with strong base.



- ❖ Boric acid (H_3BO_3) will be liberated, it is very weak acid.

- ❖ Boil the solution to remove CO_2 .


- ❖ Add mannitol because boric acid is weak acid to convert it into strong complex acid which is strong enough to be titrated directly with NaOH .

- ❖ Mannitol should be used in high conc. to prevent the hydrolysis of the complex.



sod. metaborate

❖ Mannitol should be used in high conc. to prevent the hydrolysis of this complex.



mannityl boric acid

- **Chemical factor:** *home work*
 - For borax with 0.5N HCl (0.09535 gm)
 - For borax with 1N NaOH (0.09535 gm)
 - For sod. carbonate with 0.5N HCl (0.0265gm)
 - For boric acid with 1N NaOH (0.06184gm)

Mwt. Boric acid=61.84 gm/mol

Mwt. Na₂CO₃= 106gm/mol

- **Calculations:**

we must first correct the volume of HCl & NaOH

1. If we have pure borax, the vol. of $N/2$ HCl must be equal to the vol. of $1N$ NaOH.

$V'_{\text{of } N/2 \text{ HCl}} = V'_{\text{of } 1N \text{ NaOH}}$ (the differences up to 0.3 mL)

take the average of the two volumes & calculate the % w/v of Borax.

2. Impure borax may contain Na_2CO_3 or boric acid as impurities.

a) Mixture of **Borax & Na_2CO_3**

Na_2CO_3 is basic; it can react with HCl, so both borax and carbonate will react with HCl



While NaOH react only with boric acid liberated from borax. Then after vol. correction;

$$V_1 \text{ of } N/2 \text{ HCl} > V_2 \text{ of } 1N \text{ NaOH}$$

$$V_1 - V_2 = \text{vol. of } N/2 \text{ HCl react with } \text{Na}_2\text{CO}_3$$

Calculate the % w/v for Borax and Na_2CO_3 .

b) Mixture of **Borax & boric acid**



liberated boric acid



Liberated & free boric acid

So both the liberated and free boric acid react with NaOH, then

Correct the volumes and compare:

V_2 of 1N NaOH $>$ V_1 of N/2 HCl

$V_2 - V_1 =$ vol. of 1N NaOH react with free boric acid

Calculate the % w/v for Borax and free boric acid.

Assay of Chlorinated lime

Or bleaching powder

Or Ca. Chlorohypochlorite



Assay of Chlorinated lime



- Contain not less than 30%w/w available chlorine.
- It is dull-white powder with characteristic odor, partially soluble in water and in 96% ethanol.
- Used as disinfectant.

Procedure:

The flask used in the titration is *iodine flask*

≈0.1N $\text{Na}_2\text{S}_2\text{O}_3$ →

The buret contains the titrant.

10ml of unknown
1.5gm KI
5ml of 33%v/v HAC


This flask contains the solution to be titrated and the indicator.

note: starch is added near the end point when the color of the solution changed from brown to yellow.

* **End point is colorless**



- Chemical principle:



available chlorine

The *available chlorine* is the Cl_2 which is liberated from the substance by the action of acid.

This liberated chlorine will displace an equivalent amount of I_2 from KI .



The liberated I₂ is titrated against sod. thiosulfate using starch mucilage as indicator.



Chemical factor:

1ml of 0.1N Na₂S₂O₃ \equiv 0.003545gm *available chlorine*.

Calculation:

Correct the sod. thiosulfate volume = V'

V' X ch. factor = gm of *available chlorine* in 10ml.

(Gm/10) X 100 = %w/v of *available chlorine*.

Assay of copper sulfate

Cupric sulfate $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$

Mwt. = 249.68 gm


- **Properties:** blue, fine prism or blue crystalline powder.
- **Solubility:** soluble at 20° c in 3parts of water and 3parts of glycerol, almost insoluble in 95% alcohol.
It contains not less than 98.5% and not more than 101%w/w $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$.
- **Uses:** as fungicide.

- **Procedure:** the titration should be carried in stoppered flask.

During titration the color of the solution changed from **brown** to **faint yellow** then add 1ml starch and 1gm KSCN and complete titration.

$\approx 0.1N$ **Na₂S₂O₃** →

10ml unknown
30ml D.w.
1.5 gm KI
5ml acetic acid



- **Chemical principle:**

- ❖ CuSO_4 is treated with excess KI in acidic media.

oxidizing agent

acidic media



- ❖ The formed Cupric iodide (CuI_2) is unstable.
- ❖ CuI_2 will break into cuprous iodide (CuI or Cu_2I_2) and I_2 .



unstable


yellow ppt.

- ❖ The liberated I_2 is titrated with std. $\text{Na}_2\text{S}_2\text{O}_3$ sol.



• Notes:

1. KI is added to liberate I₂.
2. Excess KI is added to dissolve I₂.
3. HAC is added to make the medium slightly acidic.
4. In acidic media the liberated I₂ is adsorbed by the ppted CuI.
5. KSCN (pot. thiocyanate) is added near the end point.



this KI will return to
react with CuSO₄ → → I₂

- **Chemical factor:**

1ml of 0.1N $\text{Na}_2\text{S}_2\text{O}_3 \equiv 0.02497\text{gm CuSO}_4 \cdot 5\text{H}_2\text{O}$

- **Calculation:**

- ❑ correct the volume of $\text{Na}_2\text{S}_2\text{O}_3 = V'$
- ❑ $V' \times \text{Ch. Factor} = \text{gm wt. of CuSO}_4 \cdot 5\text{H}_2\text{O in 10ml}$
- ❑ $(\text{Gm}/10) \times 100 = \%w/v \text{ of CuSO}_4 \cdot 5\text{H}_2\text{O in the unknown}$

Assay of Ferrous Sulfate



$M_{\text{wt.}} = 278 \text{ gm/mole}$


- **Appearance:** The blue-green heptahydrate is the most common form of ferrous sulfate ,it is bluish green crystals, freely soluble in water.
- **Content:** Ferrous Sulfate (crystal) contains 98.0—104.0% of ferrous sulfate heptahydrate ($\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$)
- **Used** in the prevention and treatment of iron deficiency.

- **Procedure:**

H₂SO₄ should be added to maintain a constant H ions concentration as explained

$\approx 0.1\text{N KMnO}_4$ 

Unknown powder
20ml of diluted H₂SO₄





The buret contains the titrant.

This flask contains the solution to be titrated and the indicator.

End point: colorless 

- **Chemical principle:**



- Ferrous sulfate is a *reducing agent*.
- It will be titrated against standard solution of potassium permanganate (*oxidizing agent*) in acidic media.
- **Note:** you must write the oxidation and reduction equations by each.
- **Chemical factor:**

Derive the ch. Factor from the previous chemical equation.

1ml. of 0.1N KMnO₄ ≡ 0.0278gm of FeSO₄

- **Calculation:**

Correct the volume of KMnO_4 consumed = V'

Then multiply the corrected volume (V') by the Ch. Factor.

to get the weight;

$V' \times \text{ch. Factor} = \text{gm of FeSO}_4 \text{ in the unknown.}$


Assay of Hydrogen Peroxide




H₂O₂: Mwt.=34.02 gm/mole

*Different H₂O₂ solution: 30%, 27%w/w and
6%, 3%w/v*


- H₂O₂ is slightly more viscous than water.
- Aqueous solution of H₂O₂ colorless and odorless. Rapidly decomposes when contact with organic matter and certain metals; Cu, Fe, Mn and if allowed to stay in alkaline solution.



- In acidic solution Fe^{2+} is oxidized to Fe^{3+} (H_2O_2 acting as an oxidizing agent):



- hydrogen peroxide acts as a reducing agent in its reaction with KMnO_4 *in our experiment*.




- **Uses:**

It is a ***disinfectant*** or ***antiseptic*** for treating wounds for its antimicrobial activity.

it is an effective ***cleaning agent*** prevents infection transmission in the hospital environment.

Procedure:

$\approx 0.1\text{N KMnO}_4$ 

10ml of unknown
25ml of D.W.
5ml of 50%v/v H₂SO₄ 



The buret contains the titrant.

This flask contains the solution to be titrated and the indicator.

End point: colorless  pink

- **Notes:**

1. We make the media acidic to;

- Prevent decomposition of H_2O_2

- Prevent reaction of $\text{KMnO}_4 \longrightarrow \text{MnO}_2$

- Increase oxidizing power of KMnO_4

2. Titration should be done slowly


- **Chemical principle:**



**You must write the oxidation and reduction equations to calculate the eqwt. for both KMnO_4 & H_2O_2 .*

- **Chemical factor:**

According to the chemical equation;



then complete the derivation of the Ch. Factor.

The Chemical factor is equal to:



- **Calculation:**
- Correct the **volume** of KMnO_4
 $N \times V = N' \times V'$ (corrected)
- Multiply the corrected volume by the Ch. Factor.
to get the **weight (wt.)** of H_2O_2
- the **%w/v** of $\text{H}_2\text{O}_2 = \text{wt.} / 10 \times 100$

- There is 2 methods to calculate the concentration of hydrogen peroxide in commercial products:

1. %W/V:

10ml of concentrated sol. was diluted to 200ml with D.W. \longrightarrow Then 20ml is taken from the dil. one and titrated with approximately 0.1N KMnO_4

<u>conc.</u>	<u>diluted</u>	$x = 20 \times 10 / 200 = \underline{1\text{ml}}$ of conc. sol. (original sol.)
10	200	
x	20	

Correct the volume of $\text{KMnO}_4 = \text{Vol.}$

Vol. x ch. Factor = wt. of H_2O_2 in 1ml.

Then calculate %w/v


2. Volume strength of the solution:

Decomposition occurs according to this equation:



Example: if the %w/v of H_2O_2 is 8.5%w/v *i.e.* It contains 8.5 gm of H_2O_2 in 100ml solution.

<u>wt.</u>	<u>Vol.</u>	
8.5	100ml	} X = 0.085 gm of H_2O_2 in 1ml
X	1ml (from the definition)	



So the volume strength of the sample is 28.

Assay of sodium
hydroxide
solution

NaOH solution

- From B.p, NaOH solution contains not less than 97.5% w/w of total alkali (as NaOH) and not more than 2.5 % w/w Na₂CO₃.
- **Assay:**
 - 10 ml of unknown(bulb pipette).
 - 25 ml of distilled water.
 - Add 2.5 ml of barium chloride solution.
 - Titrate with 1N HCl solution using 1-2 drops of phenolphthalein ind..
 - The first end point from pink → colorless(turbid)
 - To the turbid sol. add 5 drops of Bromophenol Blue ind. and complete titration with 1N HCl .
 - The second end point bluish violet → yellowish green

1N HCl solution



10 ml of unknown
25ml of D.W
2.5 ml of BaCl₂ solution
1-2 drops of ph.ph. ind. (1st titration)
5-6 drops of B.p.B ind. (2nd titration)



The buret contains
the titrant.

This flask contains the
solution to be titrated
and the indicator.

- **1st titration:**



Why HCl do not react with BaCO₃?

Why the end point is turbid?

- **2nd titration:**



definition of chemical factor: the weight of substance that is chemically equivalent to 1ml of std. solution.

- **Calculation of the *chemical factor*:**

a) From reaction of HCl with NaOH:

1Mwt of NaOH \equiv 1 Mwt HCl

1 Mwt of NaOH \equiv 1 eqwt HCL

1 \times 40 gm of NaOH \equiv 1 liter of 1N HCl

40/1000 gm NaOH \equiv 1ml of 1N HCl

0.04 gm of NaOH \equiv 1ml of 1N HCl of total
alkalinity calculated as NaOH(**chemical
factor**)

b) From reaction of 2HCl with Na_2CO_3

$2\text{Mwt of HCl} \equiv 1\text{Mwt of BaCO}_3 \equiv 1\text{Mwt Na}_2\text{CO}_3$

$1\text{Mwt Na}_2\text{CO}_3 \equiv 2\text{Mwt of HCl}$

$1\text{Mwt Na}_2\text{CO}_3 \equiv 2 \text{ eqwt of HCl}$

$\frac{1}{2} \text{ Mwt Na}_2\text{CO}_3 \equiv 1 \text{ eqwt of HCl}$

$106/2 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{liter of } 1N \text{ HCL}$

$53 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{liter of } 1N \text{ HCL}$

$53/1000 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{ml of } 1N \text{ HCL}$

$0.053 \text{ gm Na}_2\text{CO}_3 \equiv 1\text{ml of } 1N \text{ HCL (chemical factor)}$

- **Calculations :**

V_1 is the of HCl consumed in the 1st titration

V_2 is the of HCl consumed in the 2nd titration

$V_1 + V_2 = V_3$ total HCl consumed.

correct the V_3 according to this equation:

$$V \times N = V' \times N'$$

Corrected $V_3 \times 0.04 =$ gm wt. of total alkali

Then the % w/v of total alkali

Corrected $V_2 \times 0.053 =$ gm wt. of Na_2CO_3

Then the % w/v of Na_2CO_3 in the unknown

**Preparation & Standardization of ≈ 0.1 N
Potassium Permanganate (KMnO_4) solution**

Redox Reaction

◦ Preparation of $\approx 0.1N$ $KMnO_4$ solution:

- Dissolve 3.3 gm of solid $KMnO_4$ in 1000ml. distilled water, in a flask.
- Heat the solution on a hotplate for 15min. at 70-80 °c.
- Cool the solution to room temp. and allow stand for 2 days in a dark place.
- Filter the solution through the glass wool or asbestos before use.

◦ **Standardization of KMnO_4 solution:**

- ***Filter the prepared soln. through glass wool as said before.***
- ***Rinse and Fill the burette with filtered KMnO_4 soln.***
- ***Weigh accurately about 200mg sod. oxalate ($\text{Na}_2\text{C}_2\text{O}_4$), previously dried at 110°C to constant weight.***
- ***Dissolve it with 100ml of D.W. in a conical flask. Add 7ml. of conc. H_2SO_4 .***
- ***Heat the soln. to about 70°C .***
- ***Titrate against KMnO_4 soln. until a permanent pale pink color appear (persist for 15 sec.).***
- ***Record the volume of KMnO_4 consumed.***

- **OXIDATION- REDUCTION REACTIONS:**

- ❖ Methods of analysis involve a change in valence of the reacting substances, thus there is always a transfer of electrons.
- ❖ Electrons lost by reducing agent are gained by oxidizing agent. The oxidizing agent will be reduced; gain of electrons and then decreasing in oxidation number, but the reducing agent will be oxidized; losing of electrons and then increasing in oxidation number.

- We can calculate the eqwt. of the reactants:

Mwt.

$$\text{Eq. wt} = \frac{\text{Mwt.}}{\text{no. of e- gained or lost by each molecule}}$$

no. of e- gained or lost by each molecule

- **Permanganate ion** (MnO_4^-) acts as an oxidizing agent in acidic media, being reduced to Mn^{++} :



- Oxalate ion ($\text{C}_2\text{O}_4^{2-}$) acts as reducing agent:



- KMnO_4 is widely used in volumetric analysis as a standard and oxidizing agent because:
 1. It is strong oxidizing agent.
 2. The intense color of KMnO_4 solution is sufficient to signal the end point in most titrations thus it is used as self indicator.
- On the other hand,
 - 1) KMnO_4 is a strong oxidizing agent, then it will oxidize any organic matter present in D.W., so heating for 15min. is to accelerate the oxidation of these organic matter.

- 2) KMnO_4 soln. is allowed to stand for 2 days to ensure the completion of the decomposition reaction.
- 3) It is filtered through asbestos to remove all traces of manganese dioxide (MnO_2) .
- 4) KMnO_4 is unstable in the presence of direct sunlight, organic matter ex. Filter paper.
- 5) H_2SO_4 is added since KMnO_4 is better oxidizing agent in acidic media, H_2SO_4 is added to keep $[\text{H}^+]$ ion conc. constant through out the titration.
- 6) Heat to about 70°C because the oxidation of sod. oxalate is rapid enough if the temp. above 60°C .

- Sod. oxalate, not oxalic acid, is used because:
 1. The salt can be obtained in very pure condition.
 2. Free from water of crystallization.
 3. It can be dried to 130 °c without decomposition.
 4. Stable to air.

The chemical equation for the std. of pot. Permanganate:




Calculations

$$\frac{\text{KMnO}_4}{N \times V} = \frac{\text{Na}_2\text{C}_2\text{O}_4}{\frac{\text{WT}}{\text{eqwt}} \times 1000}$$

eqwt

67gm/mole



The preparation and standardization of 0.1 N sodium thiosulfate



Preparation:

Dissolve about 26gm sod. thiosulfate and 200mg sod. carbonate in 1000ml of recently boiled and cooled D.W.

Standardization:

Put into the stoppered flask the followings;

0.1–0.2 gm of potassium dichromate ($K_2Cr_2O_7$)

100ml of D.W.

3gm potassium iodide (KI)

2gm sod. Bicarbonate($NaHCO_3$)


5ml conc. HCl

- ▶ Close the flask with the stopper and swirl to mix.
- ▶ Then keep the flask in dark place for 10min.

- ▶ Rinse the stopper and inner walls of the flask with D.W.
- ▶ Titrate the liberated I_2 with sod. thiosulfate the color change from brown to faint yellow.
- ▶ Add 1 ml of starch indicator the color of the solution change to dark blue.
- ▶ Complete titration until the color change to green-blue.


≈ 0.1 N $Na_2S_2O_3$

End point: dark blue → green-blue



Notes for preparation:

- ▶ *Why do we use recently boiled and cooled water?*
 - a. to destroy sulfur bacteria
 - b. to expel CO₂



- ▶ *Na₂CO₃ is added to the solution as preservative*

Chemical principle of standardization:

- ▶ Pot. Dichromate is an oxidizing agent



- ▶ The liberated I_2 is titrated against $\text{Na}_2\text{S}_2\text{O}_3$



sod. tetrathionate

- ▶ This indirect procedure is known as *iodometry*.

Notes for standardization:

1. We use stoppered flask to prevent the loss of I_2 .
2. KI is added to liberate I_2 .
3. Excess KI is added, Why?



4. The standardization should be carried in acidic media and not in alkaline media:

because alkaline media I_2 will react with OH^- to give:



6. We put the flask in dark place for 10min.?
7. Starch ind. is added at the end (when the color of the solution change to yellowish green) or we can said the ind. is added when the iodine concentration is low because β -amylase which the soluble form of starch react :
$$\beta\text{-amylase} + \text{I}_2 + \text{I}^- \longrightarrow \text{blue colour complex}$$

the color change is reversible, the color being discharged (change to colorless) when iodine is reduced by $\text{Na}_2\text{S}_2\text{O}_3$.

This reversibility is reduced when the iodine conc. is high for this reason starch ind. should be added until most of I_2 has been reduced by $\text{Na}_2\text{S}_2\text{O}_3$.

8. The end point is a pale green and not colorless why?

Calculation:

$$\begin{array}{l} \text{Na}_2\text{S}_2\text{O}_3 \\ N \times V \\ \text{eqwt} = 49 \end{array} = \begin{array}{l} \text{K}_2\text{Cr}_2\text{O}_7 \\ (\text{Wt}/\text{eqwt}) \times 1000 \end{array}$$

The preparation and assay of Aqueous Iodine Solution (Lugol's Solution)

- ▶ It contains 5% w/v *Iodine* (limits 4.75– 5.25% w/v) and 10% w/v *Pot. Iodide* (limits 9.5– 10.5% w/v).
- ▶ Used in pre-operative treatment of thyrotoxicosis.
- ▶ Preparation:

I₂ 50gm

KI 100gm

purified water sufficient to 1000ml

- *Dissolve KI in a little amount of water then add the I₂ and shake well until the I₂ is completely dissolved, then add sufficient water to 1000ml using volumetric flask.

▶ Assay: Dilute 25ml of the original sol. to 100ml with distilled water, then: `

a) Assay of I₂,


note: we must use stoppered flask

Each 1 ml of N/10 Na₂S₂O₃ ≡ 0.01269 gm I₂

≈ N/10 Na₂S₂O₃ →

dark brown, titrate → yellow then add 1 ml starch
the color change to dark blue complete tit.
until the end point colorless

10ml of unknown
10ml D.W




b) Assay of KI,

* *titrate with M/20 KIO₃ (Pot. Iodate) until dark brown sol. change to pale brown, add 1 ml of amaranth sol.; the color of sol. Change to red and then continue tit. Until the color of sol. change to yellow -orange*

M/20 KIO₃ →

5ml of unknown
10ml D.W
20ml conc. HCl



- ▶ From the no. of ml. of M/20 KIO_3 required subtract $\frac{1}{4}$ no. of ml of N/10 $\text{Na}_2\text{S}_2\text{O}_3$ required in the assay of I_2
each 1 ml of the remainder is \equiv 0.0166 gm of KI

▶ Chemical principle for I₂ assay:

I₂ will oxidize chemically equivalent amount of sod. thiosulfate.



▶ Chemical principle for KI assay:

KIO₃ (pot. Iodate) is fairly strong oxidizing agent, it can quantitatively react with both iodide I⁻ and iodine I₂.

*If the conc. of HCl not exceed 1 N the reaction:



*In the presence of conc. HCl exceeding 4N, I₂ produced by the previous reaction is further oxidized to iodine monochloride (ICl):



* so the overall reaction between KIO₃ & KI in the presence of high conc. HCl >4N can be expressed by:



The presence of high conc. HCl is required to decrease :



stable cpd

- ▶ In this experiment we use amaranth indicator, it is a dye, dyes which are destroyed or changed in the presence of oxidizing titrant can be use as indicator.
- ▶ We can not use starch as indicator because starch_ iodide complex can not be formed in high acid conc. used in the iodate titration.

▶ ***Chemical factor. Homework***

Each 1 ml of 0.1 N $\text{Na}_2\text{S}_2\text{O}_3 \equiv 0.01269\text{gm}$ of Iodine

Each 1 ml of M/20 $\text{KIO}_3 \equiv 0.0166\text{gm}$ of KI

Each 1 ml of M/20 $\text{KIO}_3 \equiv 0.02538\text{gm}$ of iodine

▶ **Calculation:**


▶ **For I₂:**

<u>original sol.</u>	<u>Vol.</u>	}	<i>X = 2.5ml of original sol.</i>
25	100		
X	10		

V' x ch.fact. $0.01269 \times 100/2.5 = \%w/v$ of iodine in orig. sol.
V' the corrected volume of sod. thiosulfate.

▶ **For KI:**

<u>original sol.</u>	<u>Vol.</u>	}	<i>X = 1.25ml of original sol.</i>
25	100		
X	5		



Vol. of KIO₃ (req. for KI) x 0.0166 x 100/1.25 = %w/v of KI

- ▶ It is preferred to use iodate (KIO_3) of known M rather than of known N, it depends on the nature of reaction.
 - ▶ KIO_3 sol. used in the titration of sol. containing both iodide and iodine.
 - ▶ The equivalent of KIO_3 in its reaction with KI differs from its equivalent when react with iodine.
- 