

2015 - 2016



Preparation & Standardization of 0.1 N HCl Solution

Assistant Lecturer

Sahar Mohammed Shakir

Assistant Lecturer

Rana Adel Kamoon

Assistant Lecturer

Mohammed Abdul-Amir



A reagent bottle of concentrated HCl has the following informations on it's lable:

Molar mass = 36.5 g.mol⁻¹
sp.gr = 1.18
37% HCl (w/w)



How could you prepare:
1L of 0.1 N HCl solution
from this conc. HCl ?



1- Calculation of the Normality of the concentrated HCl:



$$N_{\text{HCl}} = \frac{\text{Specific gravity} * \% (w/w) * 1000}{\text{Eq.mass of HCl}}$$


$$N_{\text{HCl}} = 11.961 \text{ N}$$

2- Calculation of the volume of 11.961N HCl that should be taken to prepare 1L of 0.1N HCl soln.

$$N_1 * V_1 \text{ concentrated} = N_2 * V_2 \text{ diluted}$$

$$11.961 * V_1 = 0.1 * 1000$$

$V_1 = 8.36$ ml of concentrated HCl should be taken and diluted to the mark with distilled water in a 1000 ml volumetric flask.



Standardization of the prepared HCl solution

If the chemical is available in a pure state, e.g. anhydrous Na_2CO_3 , weigh out an exact quantity, dissolve it in water up to volume.




How could you prepare 0.5 L of 0.1 N Na_2CO_3 ? Knowing that, atomic masses of Na = 23, O = 16 and C = 12.


Substances which are not usually obtained in a pure state, e.g. mineral acids and caustic alkali, are prepared as approximate solutions and standardized against a known pure std., e.g. Na_2CO_3 as a primary std. soln.

Procedure:

- 1- Fill the burette with the prepared HCl soln.
- 2- Transfer 10 ml of exactly 0.1 N Na_2CO_3 solution (1^o - standard) in to a conical flask by using a 10 ml bulb pipette.
- 3- Add 2 drops of methyl orange as indicator. Yellow color is obtained.



4- Titrate with HCl soln. drop by drop from the burette in to the conical flask until a faint orange color is obtained.



5- The exact normality can be calculated from the following equation,




$$N_1 * V_1_{\text{Na}_2\text{CO}_3} = N_2 * V_2_{\text{HCl}}$$

Directions for reading a burette:

1- Hold an opaque card or a piece of paper behind the graduations.

2- Avoiding a parallax .
In reading volumes, the eye must be at the level of the liquid surface to avoid an error due to parallax.




3- It is common practice to use the bottom of the meniscus as the point of reference in calibrating and using volumetric equipment.

What is the reading of this burette ?

- A. 21.20 ml**
- B. 21.30 ml**
- C. 22.60 ml**
- D. 22.70 ml**



Encircle the letter that represents the right position in reading a burette .



Post Lab Exercise:

A bottle of concentrated HCl has the following informations on it's label: molar mass is 36.5 g/mol, sp.gr. 1.18 and 40% HCl (*w/w*) .

- a- What is the normality of the HCl in the bottle?***
- b- How could you prepare 2 liters of about 0.1 N HCl solution from the concentrated reagent?***