


Preparation and Standardization of 1 N NaOH Solution

- ❑ Sodium hydroxide is a strong base that is usually used to prepare standard alkaline solutions useful for volumetric analysis of acidic compounds.
- ❑ Sodium hydroxide is hygroscopic and can react with atmospheric carbon dioxide.




contaminant
(water soluble)

Preparation and Standardization of 1 N NaOH Solution

☐ preparation of 100 mL of 1 N NaOH solution

Dissolve 4.5 g of sodium hydroxide in 100 mL distilled water, allow to cool, and then add saturated barium hydroxide solution drop wise with stirring until a precipitate is formed. Leave aside allowing for complete precipitation, filter, and collect the filtrate to be standardized against 1 N HCl solution.

Preparation and Standardization of 1 N NaOH Solution



contaminant
(water soluble)



water insoluble

Preparation and Standardization of 1 N NaOH Solution

□ standardization



➤ 1 N HCl solution is used as a secondary standard

➤ phenolphthalein is used as the indicator

colourless \longrightarrow *pink*

pH: 8.3

10

Preparation and Standardization of 1 N NaOH Solution

□ procedure

- wash the burette with the D. W. and the titrant (NaOH)
- fill the burette with NaOH to a level (adjust it)
- wash a 20 mL – bulb pipette with D. W. then by a little of HCl solution; fill it to the mark with the acid
- transfer the acid into a clean conical flask; add D.W. (50 mL)
- add 2 drops of phenolphthalein indicator
- start titration by adding NaOH solution drop wise with continuous stirring until the solution changes from colourless to pink
- record the volume of NaOH solution used and calculate the normality

Note

wash the burette with water thoroughly


Preparation and Standardization of 1 N NaOH Solution

□ procedure



Preparation and Standardization of 1 N NaOH Solution

□ titration apparatus




Preparation and Standardization of 1 N NaOH Solution



**end point
(pink)**

Preparation and Standardization of 1 N NaOH Solution

□ calculations



N_1 :the normality of NaOH solution

V_1 :the volume of NaOH solution used

N_2 :the normality of HCl

V_2 : volume of HCl solution used (20mL in our experiment)

Preparation and Standardization of 1 *N* NaOH Solution

☐ Home work

Why have you used 4.5 g of NaOH to prepare 100 mL of 1 *N* NaOH solution?